


Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 1 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

LabCyte CORNEA-MODEL24
EYE IRRITATION TEST
OPERATION PROTOCOL
Ver. 2.5.6




Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 2 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

TABLE OF CONTENTS

1. RATIONALE AND BACKGROUND	4
1.1 The LabCyte CORNEA-MODEL24 EYE IRRITATION TEST	4
1.2 BACKGROUND OF LabCyte24 EIT	4
1.3 BASIS OF THE TEST METHOD	5
1.4 LIMITATION OF THE METHOD	5
1.5 BRIEF BASIC PROCEDURE	5
1.5.1 LabCyte CORNEA-MODEL24.....	6
1.5.1.1 QUALITY CONTROL OF THE TEST SYSTEMS	6
1.5.1.2 PRECAUTIONS	7
1.5.2 ASSAY QUALITY CONTROL.....	7
1.5.2.1 ASSAY ACCEPTANCE CRITERION 1: NEGATIVE CONTROL.....	7
1.5.2.2 ASSAY ACCEPTANCE CRITERION 2: POSITIVE CONTROL.....	7
1.5.2.3 ASSAY ACCEPTANCE CRITERION 3: STANDARD DEVIATION (SD).....	8
1.6 DATA INTERPRETATION PROCEDURE (PREDICTION MODEL).....	8
2. MATERIALS.....	9
2.1 LabCyte CORNEA-MODEL24.....	9
2.1.1 LabCyte CORNEA-MODEL24 KIT COMPONENTS	9
2.1.2 SHIPMENT OF LabCyte CORNEA-MODEL24	9
2.1.3 INSTRUCTIONS FOR USE OF LabCyte CORNEA-MODEL24	10
2.2 CONSUMABLES.....	10
2.3 OTHERS	10
2.3.1 EQUIPMENT/INSTRUMENTS	11
3. TEST METHOD	12
3.1 PREPARATIONS.....	12
3.1.1 POSITIVE CONTROL	12
3.1.2 NEGATIVE CONTROL.....	12
3.1.3 POLY WASH BOTTLE FOR PBS.....	12
3.2 TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT	13
3.2.1 DETECTION OF THE CHEMICALS THAT STAIN THE TISSUE	13
3.2.1.1 STEP 1 (PRELIMINARY TEST)	13
3.2.1.2 STEP 2 (FUNCTIONAL CHECK ON VIABLE TISSUE)	13
3.2.2 DETECTION OF CHEMICALS THAT DIRECTLY REDUCE WST-8	14
3.2.2.1 STEP 3 (PRELIMINARY TEST)	14
3.2.2.2 STEP 2 (FUNCTIONAL CHECK ON FREEZE-KILLED TISSUE)	15
3.3 EXECUTION OF THE TEST.....	16
3.3.1 PREPARATION OF LabCyte CORNEA-MODEL24 (DAY -1).....	16
3.3.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING (DAY 0~1)	17
3.3.2.1 PREPARATION OF WELLS FOR POST-INCUBATION (3 RD ROW)	17
3.3.2.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING.....	18
3.3.2.3 POST-EXPOSURE INCUBATION.....	20
3.3.3 APPLICATION OF SOLID TEST CHEMICALS AND RINSING (DAY 0~1).....	21
3.3.3.1 PREPARATION OF WELLS FOR SOLID TEST CHEMICAL APPLICATION (2 ND ROW) AND FOR POST-INCUBATION (3 RD ROW).....	21
3.3.3.2 APPLICATION OF SOLID TEST CHEMICALS	21

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 3 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

3.3.3.3	RINSING OF SOLID TEST CHEMICALS.....	22
3.3.4	WST-8 ASSAY (DAY 1)	24
3.3.4.1	PREPARATION OF WELLS FOR WST-8 ASSAY	24
3.3.4.2	WST-8 REACTION.....	24
3.3.4.3	SAMPLING THE REACTED WST-8 MEDIUM	25
3.3.4.4	OPTICAL DENSITY MEASUREMENTS OF THE REACTION MEDIUM	26
4.	ASSESSMENT	28
4.1	CONDITIONS FOR A VALID TEST.....	28
4.2	ASSAY CRITERIA	28
5.	References.....	29
MDS 1-1:	RECEIPT OF LabCyte CORNEA-MODEL24 (2.1.2).....	30
MDS 1-2:	TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT STEP-1 (3.2.1.1).....	31
MDS 1-3:	TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT STEP-1 (3.2.2.1).....	32
MDS 1-4:	PREPARATION OF FREEZE KILLED TISSUE (3.2.2)	33
MDS 2:	PREPARATION OF LabCyte CORNEA-MODEL24 (3.3.1).....	34
MDS 3-1(LIQUID):	APPLICATION OF LIQUID TEST CHEMICALS, RINSING AND POST-INCUBATION (3.3.2).....	35
MDS 3-2(SOLID):	APPLICATION OF SOLID TEST CHEMICALS AND RINSING (3.3.3)	37
MDS 4-1:	WST-8 ASSAY (3.3.4.1, 3.3.4.2)	39
MDS 4-2:	SAMPLING OF REACTING WST-8 SOLUTION (3.3.4.3) AND MEASUREMENT (3.3.4.4).....	40
	REVISION HISTORY	41

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 4 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

1 **1. RATIONALE AND BACKGROUND**


2 **1.1 The LabCyte CORNEA-MODEL24 EYE IRRITATION TEST**

3 The LabCyte CORNEA-MODEL24 eye irritation test (LabCyte24 EIT) is designed to identify test
4 chemicals that cause acute eye irritation by measuring cytotoxic effects using the
5 2-(2-methoxy-4-nitrophenyl)-3-(4-nitrophenyl)-5-(2,4-disulfophenyl)-2H-tetrazolium, monosodium salt
6 (WST-8) assay on the reconstructed human corneal epithelial (RhCE) model. The LabCyte24 EIT is not a
7 kit, but LabCyte CORNEA-MODEL24 tissues are commercially available at a minimum of 24 LabCyte
8 CORNEA-MODEL24 tissues per order.

9
10 **1.2 BACKGROUND OF LabCyte24 EIT**

11 Assessment of ocular irritation is an essential part of early testing procedures for the evaluation and
12 hazard classification of substances. Therefore, it plays an important role in the safety evaluation of
13 general consumer products and materials. Novel substances used in consumer products must undergo
14 comprehensive toxicological evaluation for eye irritation and a variety of other adverse outcomes. To
15 date, the *in vivo* Draize eye test has been the international standard assay for acute ocular toxicity
16 evaluation (irritation and corrosion) and is described (including optimizations and refinements) in OECD
17 Test Guideline 405. However, the use of this test has been questioned and strongly criticized for ethical
18 concerns related to animal welfare, because it is painful to the rabbits. Thus, alternative strategies and
19 tests are urgently required in order to evaluate the eye irritation potential of new chemicals. Corneal
20 epithelial cells on the surface of the eye are first to be exposed to substances and have been widely
21 studied for links to the biological role of tissue and gene regulation. Three-dimensional RhCE models are
22 useful as a multilayered, standardized tissue that mimics the human corneal epithelium.

23 The LabCyte24 EIT was developed as a replacement for the Draize eye irritation test. The Draize
24 scoring system is heavily weighted towards corneal damage (80 out of 110 total score), because
25 irreversible damage to the cornea can lead to blindness. Since damage to the cornea is so crucial for
26 human health, corneal tissue can be considered a useful tool for the development of *in vitro* eye irritation
27 testing.

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 5 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

29 **1.3 BASIS OF THE TEST METHOD**

30 Damage induced by eye irritants generally progresses from the corneal epithelium through the stroma
31 and potentially to the endothelium, and the LabCyte24 EIT is able to provide information on the first
32 stages of this progression. Irritants damage cells while penetrating the corneal epithelium layer, and the
33 cytotoxic progression can be estimated by analyzing cell viability in the LabCyte CORNEA-MODEL24
34 tissue using standardized methods. Although the tissues represent only the corneal epithelium, (very mild
35 responses would also be reflective of some conjunctival irritation), it can be used to estimate deeper
36 damage as far as the stroma by the analyzing cell viability.


37 The relative viability of the tissue exposed to a test chemical is measured using the WST-8 assay
38 immediately after exposure and again after a post-exposure period. A viability of 40% of the negative
39 control value was used as the cutoff in identifying test chemicals as an irritant (GHS Category 1 or 2) or
40 an non-irritant (GHS No Category). Some culture environments might allow the detection of very small
41 quantities of cytokines secreted by the corneal epithelial tissue in response to topical application of test
42 chemicals.

43
44 **1.4 LIMITATION OF THE METHOD**

45 One limitation of this assay method is potential interference of the test chemical with the WST-8
46 endpoint. A colored test chemical or one that directly reduces WST-8 (and thereby mimics dehydrogenase
47 activity of the cellular mitochondria) could interfere with the WST-8 endpoint. This is only a problem,
48 however, if there are significant residual levels of the test chemical on or in the tissue at the time of the
49 WST-8. Although this is considered unlikely, if it did happen, both the actual metabolic WST-8 reduction
50 and the quasi direct WST-8 reduction by a colored test chemical can be quantified using the procedure
51 described in **Section 3.2 “TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH**
52 **WST-8 ENDPOINT”**.

53
54 **1.5 BRIEF BASIC PROCEDURE**

55 On the day of receipt, LabCyte CORNEA-MODEL24 tissues are incubated overnight to release
56 transport–stress-related compounds and debris.

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 6 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

57 For liquid test chemicals, tissues are topically exposed for 1 minute. Preferably, three tissues are used
58 for each test chemical as well as for the positive and negative controls. After exposure, tissues are
59 thoroughly rinsed and blotted to remove the test chemical, then transferred to fresh medium and
60 post-incubated for 24 hours. For solid test chemicals, tissues are exposed for 24 hours but are not
61 post-incubated.

62 After post-incubation of tissue exposed to liquid test chemicals or after exposure of tissue exposed to
63 solid test chemicals, the tissues are each transferred to a well containing the WST-8 medium in a 1:10
64 dilution with Earle balanced salt solution (EBSS). After a four-hour WST-8 incubation, the orange
65 water-soluble formazan salt is formed in the WST-8 medium by cellular mitochondria and optical density
66 (OD) of the WST-8 medium is measured using a spectrophotometer at 450 nm and 650 nm as reference.
67 Relative cell viability is calculated for each tissue as % of the mean of the negative control tissues. Test
68 chemicals that produce a relative cell viability below 40% of the negative control are predicted to be
69 irritants.

70


71 **1.5.1 LabCyte CORNEA-MODEL24**

72 The LabCyte CORNEA-MODEL24 is a commercially available RhCE model produced by Japan
73 Tissue Engineering Co. Ltd. It comprises normal human corneal epithelial cells that are derived from a
74 human cornea. The cells are cultured with 3T3-J2 cells as a feeder layer in order to expand them while
75 maintaining their phenotype (Rheinwald and Green, 1975; Green, 1978). Reconstruction of the human
76 cultured corneal epithelial tissue is achieved by cultivating and proliferating the corneal epithelial cells
77 on an inert filter substrate with a surface area of 0.3 cm² at an air-liquid interface for 13 days using an
78 optimized medium containing 5% fetal bovine serum. This results in the formation of a multilayer
79 structure comprising a fully differentiated corneal epithelium with features mimicking those of a normal
80 human corneal epithelium. For delivery, LabCyte CORNEA-MODEL24 tissues are embedded in an
81 agarose gel containing a nutrient solution and shipped in 24-well plates.

82

83 **1.5.1.1 QUALITY CONTROL OF THE TEST SYSTEMS**

84 LabCyte CORNEA-MODEL24 tissue is manufactured in accordance with well-defined quality

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 7 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

85 assurance procedures. Each production batch was comes with quality control documentation that
86 identifies storage conditions, RhCE instructions for use, lot number and origin, histology (demonstration
87 of human multilayered corneal epithelial-like structure), cell viability, and barrier function integrity
88 ($0.1 \leq IC_{50} \leq 0.4$).

89
90

91 **1.5.1.2 PRECAUTIONS**

92 Corneal epithelial cells are taken from healthy donors who are free of HIV or hepatitis. Nevertheless,
93 always adhere to the following procedures for the handling of biological materials:

- 94 (a) Always wear gloves during handling of the eye and kit components.
95 (b) All corneal epithelial tissue as well as all material and media that came in contact with it should be
96 decontaminated prior to disposal using special containers or autoclaving.

97

98 **1.5.2 ASSAY QUALITY CONTROL**

99

100 **1.5.2.1 ASSAY ACCEPTANCE CRITERION 1: NEGATIVE CONTROL**

101 The absolute OD of the negative control for either liquid or solid test chemicals (NC-L or NC-S)
102 tissues (treated with sterile PBS for liquid test chemicals or untreated for solid test chemicals) in the
103 WST-8 assay is an indicator of tissue viability obtained in the testing laboratory after shipping and
104 storing procedures and under specific conditions of use.

105 **$0.5 \leq \text{Mean OD (A450/650) measured value} \leq 1.3$**


106

107 **1.5.2.2 ASSAY ACCEPTANCE CRITERION 2: POSITIVE CONTROL**

108 Ethanol is used as the positive control (PC) for liquid test chemicals (PC-L) and is tested concurrently
109 with the liquid test chemicals. Lauric acid is used as the PC for solid test chemicals (PC-S) and is tested
110 concurrently with the solid test chemicals.

111 Concurrent here means that the PC-L and the PC-S are to be tested for each run.

112 **Mean tissue viability $\leq 40\%$**

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 8 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

113

114 **1.5.2.3 ASSAY ACCEPTANCE CRITERION 3: STANDARD DEVIATION (SD)**

115 Since eye irritation potential is predicted from the mean viability of three individual tissues, the
116 variability of tissue replicates must kept at an acceptably low level.

117 **Standard Deviation (SD) of tissue viability of three identically treated replicates for the negative**
118 **control, positive control, and test chemicals \leq 18%**

119

120 **1.6 DATA INTERPRETATION PROCEDURE (PREDICTION MODEL)**

121 The United Nations Globally Harmonized System (GHS) (United Nations, 2003) was used as a
122 reference for the in vivo eye irritation classification of test chemicals.

123 For the purpose of this EIT, an eye irritant is defined as a substance that induces reversible ocular
124 lesions after administration to rabbits.


125 According to GHS classifications, a substance is an irritant (Category 1 or 2) if the mean relative
126 viability of three individual tissues exposed to the test chemical is falls below 40% of the mean viability
127 of the negative control. **(Refer to Table 1.)**

128 **Table 1** Prediction model of LabCyte24 EIT

<i>In vitro</i> results	Prediction
Tissue viability is \leq 40%	Category 1 or 2
Tissue viability is $>$ 40%	No Category

129

130

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 9 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

131 **2. MATERIALS**

132 **2.1 LabCyte CORNEA-MODEL24**

133 **2.1.1 LabCyte CORNEA-MODEL24 KIT COMPONENTS**

134 LabCyte CORNEA-MODEL24 kit components are shown in **Table 2**.

135 **Table 2** LabCyte CORNEA-MODEL24 Kit Components

Component	Qty	Description
LabCyte CORNEA-MODEL24 plate	1 plate	Contains 24 culture inserts with tissues fixed in nutritive agar medium for transport (usable area: 0.3 cm ²).
Assay Medium	1 bottle	Basic medium for incubation (30 mL). Keep refrigerated.
24-well plate	1 plate	Blank plate for use in assay. Store at room temperature.

136


137 **2.1.2 SHIPMENT OF LabCyte CORNEA-MODEL24**

138 LabCyte CORNEA-MODEL24 is packed in a special Icompo container available from and delivered
139 by NIPPON EXPRESS CO., LTD. After the kit is delivered, examine the contents and make sure that all
140 components (LabCyte CORNEA-MODEL24 plate, assay medium, and 24-well assay plate) are included
141 in the package. Also confirm lot numbers and expiration dates. Record details in Methods Documentation
142 Sheet (MDS) 1. (See **MSD-1**).

143 NIPPON EXPRESS will pick up the Icompo container at a later date (generally, the day after delivery),
144 and it should be returned together with an shipping invoice and the insulating materials.

145

146

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 10 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

147 **2.1.3 INSTRUCTIONS FOR USE OF LabCyte CORNEA-MODEL24**

148 Always incubate all of the culture inserts after opening the package. Do not store the culture inserts
149 after opening.

150 The human corneal epithelial tissue cells used in the LabCyte CORNEA-MODEL24 originate from a
151 normal human donor and are HIV-, HBV-, HCV-, and HPV-negative. They are to be handled,
152 nevertheless, with due care and in accordance with laboratory biosafety guidelines for handling
153 human-derived materials.

154

155 **2.2 CONSUMABLES**

156 The following consumables are required.

157 *The following quantities are necessary to assay between one and six 6 test chemicals at a time.

- 158 • Assay Medium, 100 mL (J-TEC: 402250) 1
- 159 • Cell Counting Kit-8, 500 test (Dojindo: CK04) 4
- 160 • 24-well assay plate (Becton, Dickinson and Company: 353047) 7-8
- 161 • 96-well plate (Becton, Dickinson and Company: 353072 or equivalents) 1
- 162 • PBS, 500 mL (Invitrogen: 14190-144 or equivalents) 2 or 3 bottle
- 163 • Earle balanced salt solution (EBSS), 500 mL (SIGMA-ALDRICH: E3024) 1 bottle
- 164 • Sterile cotton buds (JAPAN COTTON BUDS: 10A754D or equivalents) 1 box
- 165 • Micro-pipette tips (sterile: 10~200µL, 200-1000µL)
- 166 • Microtubes (1.5mL)
- 167 • Dish (10cm)
- 168 • Paper towel

169


170 Convenient consumable items are shown followings.

171 Also, it would be convenient to have the following.

- 172 • Capillary & piston for positive-displacement-type pipette (10-100µL)

173

174 **2.3 OTHERS**

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 11 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

175 **2.3.1 EQUIPMENT/INSTRUMENTS**


- 176 • Safety cabinet (or clean bench)
- 177 • Water bath (37°C)
- 178 • CO₂ incubator (37°C, 5% CO₂, capable of maintaining high humidity)
- 179 • Autoclave
- 180 • 96-well multi-plate reader (required filters: 450 nm, 650 nm)
- 181 • Precision balance (0.1 mg)
- 182 • Aspirator
- 183 • Stop-watches
- 184 • Adjustable micro-pipette (10–200 µL, 200–1000 µL)
- 185 • Sharp-edged forceps (sterile)
- 186 • Micro spatula (sterile)
- 187 • Beaker (1–2 L: sterile)
- 188 • Sterilizable poly wash bottle (500–1000 mL: sterile) with wide mouth (mouth > 3-mm dia.)
- 189 • Mortar with pestle

190

191 Also, it would be convenient to have the following.

- 192 • Positive-displacement-type pipette (10–100 µL)

193

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 12 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

194 **3. TEST METHOD**

195 *Procedures described in Sections 3.1.1 to 3.1.3 and Sections 3.3.1 to 3.3.3 are to be performed
196 aseptically in a safety cabinet or clean bench. Procedures other than those mentions in the previous
197 sentence need not be performed aseptically. Refer to **Section 2.1.3 “INSTRUCTIONS FOR USE**
198 **OF LabCyte CORNEA-MODEL24”**.

199
200 **3.1 PREPARATIONS**

201 **3.1.1 POSITIVE CONTROL**

- 202 (1) Use Ethanol as a positive control for liquid test chemicals.
203 (2) Use Lauric acid as a positive control for solid test chemicals.


204
205 **3.1.2 NEGATIVE CONTROL**

- 206 (1) Use PBS as a negative control for liquid test chemicals.
207 (2) Use non-treatment as a negative control for solid test chemicals.

208
209 **3.1.3 POLY WASH BOTTLE FOR PBS**

- 210 (1) Sterilize poly wash bottle using an autoclave.
211 (2) Fill the sterilized poly wash bottle with sterile PBS.

212
213
214

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 13 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

215 **3.2 TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT**

216 There are two kinds of test chemicals that interfere with the WST-8 assay.

217 (a) Test chemicals that stain corneal epithelial tissues directly.

218 (b) Test chemicals that react directly with WST-8.

219 Test chemical that stain the corneal epithelial tissues could possibility transfer from the corneal
220 epithelial tissues to the WST-8 reaction buffer and affect the OD measurements. Test chemicals that react
221 directly with WST-8 can also affect OD measurements, if the test chemical is present in the corneal
222 epithelial tissues when the WST-8 viability test is performed. The procedure for detecting such test
223 chemicals is described below.

224

225 **3.2.1 DETECTION OF THE CHEMICALS THAT STAIN THE TISSUE**

226 **3.2.1.1 STEP 1 (PRELIMINARY TEST)**

227 (1) Add 50µL (Liquid) or 10mg (Solid) of the test chemical into wells of 24-well assay plate
228 preliminarily filled with 0.5mL of distilled water. Untreated distilled water is used as control.

229 (2) Close the lid of 24-well assay plate and incubate the mixture in CO₂ incubator for 15 minutes.

230 (3) After incubation, shake the mixture gently and evaluate the staining of the distilled water
231 macroscopically.

232 (4) When the color of the solution changes significantly, the test chemical is presumed to have the
233 potential to stain the tissue and a functional check on viable tissues (Step2) should be performed.

234 When the color of the solution does not change significantly, the test chemical is determined not to
235 have a potential to stain the tissue.


236 (5) Record the details of steps 1 to 4 above in **MDS 1-2**.

237

238 **3.2.1.2 STEP 2 (FUNCTIONAL CHECK ON VIABLE TISSUE)**

239 (1) dd 50 µL (liquid) or 10 mg (solid) of the test chemical, which clearly changed the color of the
240 distilled water in step 1, onto the surface of the epidermis tissues. Distilled water is used as negative
241 control.

242 (2) Follow all procedures described in the Section **3.3 EXECUTION OF THE TEST**. However,
243 incubate the tissue for 3 hours in medium (EBSS) without WST-8 instead of incubating in the
244 diluted WST-8 medium to evaluate the staining of the WST-8 medium (correcting tissue).

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 14 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

245 (3) The corrected OD is calculated using the following formula.

246 **Corrected OD = A – (B – C),**

247 where

248 A is the OD of viable tissue exposed to a test chemical,

249 B is the mean OD of correcting tissue exposed to a test chemical, and

250 C is the mean OD of correcting tissue exposed to the negative control.

251

252 (4) If a corrected OD is below 0, the OD is considered to be 0.

253 (5) When a cell viability that is calculated according to the procedures described in **Section 3.3.4.4** is
254 $\leq 40\%$, the test chemical is predicted to be an irritant (GHS Category 1 or 2) and there is no need
255 to calculate a corrected value.

256

257

258 **3.2.2 DETECTION OF CHEMICALS THAT DIRECTLY REDUCE WST-8**

259 **3.2.2.1 STEP 3 (PRELIMINARY TEST)**

260 (1) Dilute the cell counting kit-8 (WST-8 stock solution) with EBSS (Cell Counting Kit-8:EBSS = 1:10),
261 and then prepare the diluted WST-8 medium.


262 *Prepare immediately prior to use and use as soon as possible.


263 Fill each well of the 24-well assay plate with 0.3 mL of diluted WST-8 medium.

264 (2) Add 50 μ L of a liquid test chemical or 10 mg of a test chemical to the wells of 24-well assay plate.
265 The diluted WST-8 medium is used as control.

266 (3) Close the lid of 24-well assay plate and incubate the in a CO₂ incubator for about 4 hours.

267 (4) After incubation, shake the mixture gently and evaluate the staining of the diluted WST-8 medium
268 macroscopically.



Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 15 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

269 **Photo 1** Example of test chemicals to interfere the WST-8 assay directly (STEP 1).

270 When test chemicals like m-Phenylene diamine, Monoethanolamine, or Calcium thioglycollate have
271 colored the diluted WST-8 medium, Step 2 must be performed.

272 (5) Significant coloring of the diluted WST-8 medium by the test chemical can interfere with the WST-8
273 assay and therefore requires the additional functional check of Step 2.

274 (6) Record the details of steps 1 to 5 above in **MDS 1-3**.

275

276 **3.2.2.2 STEP 4 (FUNCTIONAL CHECK ON FREEZE-KILLED TISSUE)**

277 (6) Add 50 µL of a liquid test chemical or 10 mg of solid a test chemical that clearly changed the color
278 of the diluted WST-8 medium (**3.2.2.1. STEP 3**) to the surface of the corneal epithelial tissues.

279 (7) Follow all procedures described in **Section 3.3 “EXECUTION OF THE TEST”**. Instead of using
280 viable corneal epithelial tissues, however, use corneal epithelial tissues that were freeze-killed at
281 -80°C or lower for 30 min.

282 Record the details of freeze-killing the tissue in **MDS 1-3**.

283 (8) The corrected OD is calculated using the following formula.

284 **Corrected OD = A – (B – C),**

285 where

286 A is the OD of viable tissue exposed to a test chemical,

287 B is the mean OD of freeze-killed tissue (correcting tissue) exposed to a test chemical, and

288 C is the mean OD of freeze-killed tissue (correcting tissue) exposed to the negative control.


289

290 (9) If a corrected OD is below 0, the OD is considered to be 0.

291 (10) When a cell viability that is calculated according to the procedures described in **Section 3.3.4.4** is
292 ≤ 40%, the test chemical is predicted to be an irritant (GHS Category 1 or 2) and there is no need
293 to calculate a corrected value.

294

295

Version 2.5.6	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 16 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

296 **3.3 EXECUTION OF THE TEST**

297

298 **3.3.1 PREPARATION OF LabCyte CORNEA-MODEL24 (DAY -1)**

299 (1) Warm the assay medium to 37°C for 30 minutes
300 in a water bath.

301 (2) Fill the six wells of the 1st row of each assay
302 plate for LIQUID/SOLID with 0.5 mL/well of
303 the warm assay medium.

304 **See Fig. 1.**

305 (3) Open the LabCyte CORNEA-MODEL24
306 aluminum package.

307 (4) Open the LabCyte CORNEA-MODEL24 plate lid and pick up the culture inserts using sterile
308 forceps.

309 *Do not touch the surface of the corneal epithelial tissue in the culture inserts.

310 *Use forceps to remove any agar medium
311 sticking to the outside of the culture inserts.

312 (5) Transfer the culture inserts to the assay medium
313 in the wells of the 1st row using sterile forceps.

314 **See Fig. 2.**

315 *Do not allow air bubbles to form under the
316 tissue inserts.

317 (6) Close the lid on the plate and place it in a CO₂ incubator.



318 (7) Incubate overnight (15–30 hours) until ready to perform **Section 3.3.2 “APPLICATION OF**
319 **LIQUID TEST CHEMICALS, RINSING AND POST-INCUBATION”** or **Section 3.3.3**
320 **“APPLICATION OF SOLID TEST CHEMICALS AND RINSING”**.


321 (8) Record the details of steps 1 to 7 above in MDS 2.

322

323

324



Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 17 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

325 **3.3.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING (DAY 0~1)**


326

327 **3.3.2.1 PREPARATION OF WELLS FOR POST-INCUBATION (3RD ROW)**

328 (1) Warm the assay medium to 37°C for 30 minutes
329 using a water bath.

330 (2) Remove the assay plate for LIQUID from the CO₂
331 incubator.

332 (3) Open the lid of the assay plate for LIQUID, and
333 use a micropipette to fill the six wells of the 3rd
334 row with 0.5 mL/well of the warm assay medium.



335 **See Fig. 3.**


336 (4) Close the lid of the assay plate and perform **Section 3.3.2.2 “APPLICATION OF LIQUID TEST**
337 **CHEMICALS AND RINSING”** continuously.

338 (5) If the test chemicals are not applied immediately, store the Assay Plate for LIQUID in a CO₂
339 incubator until ready apply but for no more than 12 hours.

340 (6) Record the details of steps 1 to 5 above in **MDS 3-1**.

341


342

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 18 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

343 **3.3.2.2 APPLICATION OF LIQUID TEST CHEMICALS AND RINSING**

344 (1) Remove the assay plate for LIQUID from the CO₂ incubator.

345 (2) Using a micropipette, apply 50 µL of a liquid test
346 chemical to the surface of the corneal epithelial
347 tissues in the 1st row of the assay plate. Each test
348 chemical is to be tested in three wells (N=3).
349 Carefully apply the test chemical to the central
350 part of each corneal epithelial tissue. After
351 application, close the lid of the assay plate and tap
352 the sides of the plate to spread the liquid test
353 chemicals to spread out over the entire corneal




354 epithelial surface. If necessary, use a micro spatula to spread the liquid test chemical over the entire
355 surface. Take care not to press down on the surface of the corneal epithelial with the spatula.

356 *For viscous LIQUID test chemicals, use a wide orifice cell saver tip (**See Photo 2.**) or positive
357 displacement type pipette.

358 *Use a pipette or other equipment to familiarize
359 yourself beforehand with the characteristics
360 the test chemicals.

361 *Assay no more than two test chemicals on one
362 24-well assay plate.




363 **See Fig. 4.**

364 Each chemical is tested in three wells, using
365 three tissues (N = 3).

366 (3) Apply test chemicals to each well at an interval of one to three minutes.

367 (4) Close the lid and incubate each well for 60 ±10 seconds in the cabinet at room temperature.

368 *Keep the lid of the assay plate closed at all times except when applying test chemicals. Leaving the
369 lid open could affect the quantity of the test chemical in the well due to air circulation in the
370 cabinet.

Version 2.5.6	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 19 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

371 (5) Wait for 60 ±10 seconds after applying a test chemical,
372 then open the assay plate for LIQUID and pick up a
373 culture insert with sterile forceps.

374 (6) Discard the test chemical on the tissue by tilting the
375 insert and tapping it on a beaker. Fill the culture insert
376 to overflowing with PBS from a poly wash bottle.

377 *Apply a stream of PBS directly to the tissue surface to
378 wash away the test chemical. Maintain a uniformly
379 strong stream of PBS while washing the tissue.

380 **See Photo 3.**

381 *To avoid damaging the tissue with too forceful a
382 stream, use a wide-mouth nozzle on the poly wash
383 bottle.

384 (7) Tilt the insert to discard the PBS into the beaker.
385 Remove as much of the PBS inside the culture insert as
386 possible by tapping it on the beaker.

387 **See Photo 4.**

388 (8) Repeat steps 6 and 7 at least 10 times to remove as much as possible of the residual test chemical on
389 the tissue surface.

390 *Depending upon the physical properties of the test
391 chemical, bubbles might form in an insert during
392 washing. Continue washing until all bubbles
393 disappear.

394 (9) Using a sterile cotton bud, gently remove as much as
395 possible of the leftover PBS both inside and outside the
396 culture insert.

397 **See Photo 5.**

398 *Take care not to press down on the surface of the tissue

Photo 3 Rinse 1




Photo 4 Rinse 2



Photo 5 Rinse 3



Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 20 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

399 with the sterile cotton bud.

400 (10) Repeat steps 6 to 9 to remove any residual test chemical from the corneal epithelial tissue surface.

401 (11) Blot the culture insert and then place it in the well
402 of the same column in the 3rd row.

403 **See Fig. 5.**

404 *Take care to prevent air bubbles from forming
405 under the culture inserts.

406 (12) Repeat steps 1 to 11 for all culture inserts at one- to
407 three-minute intervals.

408 (13) Record the details of steps 1 to 12 above in **MDS 3-1**.

409

410 3.3.2.3 POST-EXPOSURE INCUBATION

411 (1) After performing **Section 3.3.2.2 “APPLICATION OF LIQUID TEST CHEMICALS AND**
412 **RINSING”**, close the lid of the assay plate for LIQUID and place it in a CO₂ incubator as soon as
413 possible.


414 (2) Incubate for 24 ±1 hours.


415 (3) Record the details of steps 1 and 2 above in **MDS 3-1**.

416

417

418



Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 21 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

419 **3.3.3 APPLICATION OF SOLID TEST CHEMICALS AND RINSING (DAY 0–1)**

420


421 **3.3.3.1 PREPARATION OF WELLS FOR SOLID TEST CHEMICAL APPLICATION (2ND**
422 **ROW) AND FOR POST-INCUBATION (3RD ROW)**

423 (1) Warm the assay medium to 37°C for 30 minutes using a water bath.

424 (2) Remove the assay plate for SOLID from the CO₂
425 incubator.

426 (3) Open the lid of the assay plate for SOLID and use
427 a micropipette to fill the 12 wells in the 2nd and 3rd
428 rows with 0.5 mL/well of warm assay medium .

429 **See Fig. 6.**



430 (4) Close the lid of the assay plate and perform

431 **Section 3.3.3.2 “APPLICATION OF SOLID TEST CHEMICALS”** continuously.

432 (5) If the test chemicals are not applied immediately, store the Assay Plate for LIQUID in a CO₂
433 incubator until ready apply but for no more than 12 hours.

434 (6) Record the details of steps 1 to 5 above in **MDS 3-2**.

435

436 **3.3.3.2 APPLICATION OF SOLID TEST CHEMICALS**

437 (1) Remove the assay plate for SOLID from the CO₂
438 incubator.

439 Using a precision balance, weigh out 10 ±2 mg of
440 the solid test chemicals. If necessary, crush and
441 grind the solid test chemicals in a mortar with
442 pestle. Apply the solid test chemical to the surface
443 of the corneal epithelial tissue.

444 **See Photo 6.**

445 If necessary, use a micro spatula to spread the test chemical gently over the entire surface. Use three
446 wells per test chemical (N = 3).



447 **See Fig. 7.**
448 Each chemical is tested in three wells, using
449 three tissues (N = 3).

450 (2) Place the exposed culture insert in the well of
451 the same column in the 2nd row (chemical
452 application).

453 **See Fig. 8.**

454 *Take care to prevent air bubbles from forming
455 under the culture inserts.

456 (3) Close the lid of the assay plate for SOLID and
457 place it in a CO₂ incubator. Incubate for
458 24 ±1 hours.

459 (4) Record the details of steps 1 to 4 above in
460 **MDS 3-2.**

461

462 **3.3.3.3 RINSING OF SOLID TEST CHEMICALS**

463 (1) After incubation, remove the assay plate for SOLID from the CO₂ incubator.

464 (2) Open the assay plate for SOLID and pick up a culture insert with sterile forceps.



465 (3) Discard test chemicals on the tissue by tilting the insert and tapping it on a beaker. Fill the culture
466 insert to overflowing with PBS from a poly wash bottle.

467 *Apply a stream of PBS directly to the tissue surface to
468 wash away the test chemical. Maintain a uniformly
469 strong stream of PBS while washing the tissue.

470 **See Photo 7.**

471 *To avoid damaging the tissue with too forceful a
472 stream, use a wide-mouth nozzle on the poly wash
473 bottle.

474 (4) Tilt the insert to discard the PBS into the beaker.



461


467 *Apply a stream of PBS directly to the tissue surface to
468 wash away the test chemical. Maintain a uniformly
469 strong stream of PBS while washing the tissue.

470 **See Photo 7.**

471 *To avoid damaging the tissue with too forceful a
472 stream, use a wide-mouth nozzle on the poly wash
473 bottle.

474 (4) Tilt the insert to discard the PBS into the beaker.



Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 23 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

475 Remove as much of the PBS inside the culture insert as
476 possible by tapping it on the beaker.

477 **See Photo 8.**

478 (5) Repeat steps 3 and 4 at least 10 times to remove as
479 much as possible of the residual test chemical on the
480 tissue surface.

481 (6) Using a sterile cotton bud, gently remove as much as
482 possible of the leftover PBS both inside and outside
483 the culture insert.

484 **See Photo 9.**

485 (7) If it proves difficult to remove completely all the
486 residual test chemical from the corneal epithelial tissue
487 surface, remove as much as possible and continue to
488 step 8.

489 (8) Place the rinsed culture insert in the well of the same
490 column in the 3rd row (rinse).


491 **See Fig. 9.**

492 *Take care to prevent air bubbles from forming
493 under the culture inserts.

494 (9) Record the details of steps 1 to 8 above in **MDS**
495 **3-2.**

496 After rinsing step, perform **Section 3.3.4**
497 **“WST-8 assay”** continuously.

498
499



500 **3.3.4 WST-8 ASSAY (DAY 1)**

501

502 **3.3.4.1 PREPARATION OF WELLS FOR WST-8 ASSAY**

- 503 (1) Warm the EBSS to 37°C for 30 minutes using a water bath.
- 504 (2) Dilute the cell counting kit-8 (WST-8 stock solution) with EBSS (Cell Counting Kit-8:EBSS = 1:10),
- 505 and then warm the diluted WST-8 medium. Blanks for the WST-8 assay are prepared using the
- 506 diluted WST-8 medium as follows.

507 *Prepare immediately prior to performing the

508 WST-8 assay.

- 509 (3) Remove the assay plate for LIQUID from the CO₂
- 510 incubator or prepare the assay plate for SOLID.

- 511 (4) Open the lid of the assay plate and use a
- 512 micropipette to fill each well of the 4th row with
- 513 0.3 mL/well of the warm diluted WST-8 medium.

514 **See Fig. 10.**

515 Close the lid of the assay plate and perform **Section 3.3.4.2 “WST-8 REACTION”** continuously.



- 516 (5) Record the details of steps 1 to 4 above in **MDS 4-1**.

517

518 **3.3.4.2 WST-8 REACTION**

- 519 (1) Add about 20 mL of PBS each to two dishes: PBS dish 1 and PBS dish 2.
- 520 (2) Open an assay plate (for either LIQUID or SOLID) and pick up a culture insert with sterile forceps.
- 521 (3) Clean the bottom of the culture insert so it is free
- 522 of residual culture medium by washing it first in
- 523 PBS dish 1 and the in PBS dish 2. After washing,
- 524 wipe the bottom with a paper towel.
- 525 (4) After washing and wiping, place the culture insert
- 526 in the well of the same column in the 4th row.

527 **See Fig. 11.**



528 *Take care to prevent air bubbles from forming under the culture inserts.

529 (4) Close the lid of the assay plate and place it in the CO₂ incubator.

530 Incubate for 4 hours ± 20 minutes.

531 Record the details of steps 1 to 5 above in **MDS 4-1.**

532


533 **3.3.4.3 SAMPLING THE REACTED WST-8 MEDIUM**

534 (1) After incubation, remove the assay plate from the CO₂ incubator.

535 (2) Open the lid of the assay plate and remove the culture inserts from the 4th row with forceps.

536 Transfer 200 µL of the reacted WST-8 dilution medium into the wells of a 96-well plate.

537 *Figs. 12A and 12B show typical allocations on a 96-well plate for both living and correcting tissue.



538

539

540

541 (3) Record the details of steps 1 to 2 above in **MDS 4-2**.

542

543 **3.3.4.4 OPTICAL DENSITY MEASUREMENTS OF THE REACTION MEDIUM**

544 (1) Using a 96-well plate reader, measure OD at 450 nm and 650 nm and then used the following
545 equation to determine a composite OD for each well.

546


547
$$\text{Composite OD} = (\text{OD}_{\text{TC}} \text{ at } 450 \text{ nm} - \text{OD}_{\text{blank}} \text{ at } 450 \text{ nm}) - (\text{OD}_{\text{TC}} \text{ at } 650 \text{ nm} - \text{OD}_{\text{blank}} \text{ at } 650 \text{ nm})$$

548 *If the plate reader can be programmed to perform this calculation automatically, then only the
549 composite OD value need be recorded.

550

551 (2) Calculate the Mean OD for the negative control, a cell viability for each individual tissue, and a
552 mean cell viability (including SD) for each test chemical using the following equations.

553

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 27 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

$$\text{Mean OD}_{\text{NC}} = \frac{\text{Sum of the OD}_{\text{NC}} \text{ for three replicate tissues}}{3}$$

554

$$\text{Tissue cell viability (\%)} = \frac{\text{Each tissue OD}_{\text{TC}}}{\text{Mean OD}_{\text{NC}}} \times 100$$

555


$$\text{Mean cell viability (\%)} = \frac{\text{Sum total of cell viability (\%)} \text{ for three replicate tissues}}{3}$$

556

557 (3) Record the details of steps 1 and 2 above in **MDS 4-2**.

558

559

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 28 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

560 **4. ASSESSMENT**

561

562 **4.1 CONDITIONS FOR A VALID TEST**

563 An eye irritation test is considered valid if all three of the following criteria have been met.

564

- 565 • Tissue viability: $0.5 \leq \text{mean OD (A450/650) measured value for negative control} \leq 1.3$
- 566 • Positive control: mean tissue viability for positive control $\leq 40\%$
- 567 • SD: SD (negative control, positive control and each test chemicals) of tissue viability of 3
- 568 identically treated replicates $\leq 18\%$

569

570 **4.2 ASSAY CRITERIA**

571 The criteria for in vitro prediction are shown below.

572 After exposure to a chemical, if cell viability is 40% or less, the chemical is predicted to be an
573 irritant (GHS Category 1 or 2), otherwise it is predicted to be a non-irritant (GHS No Category).

574 See Table 3.


575

576 **Table 3** Prediction model of LabCyte24 EIT

Tissue Viability	Prediction
$\leq 40\%$	Category 1 or 2
$> 40\%$	No Category

577

578


Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 29 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

579 **5. References**

580 Green H. (1978) Cyclic AMP in relation to proliferation of the epidermal cell: New view. *Cell*, 15,
581 801-811.

582 Rheinwald J.G. and Green H. (1975) Serial cultivation of strains of human epidermal keratinocytes:
583 The formation of keratinizing colonies from single cells. *Cell*, 6, 331-343.

584

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 30 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

585 **MDS 1-1:**
586 **RECEIPT OF LabCyte CORNEA-MODEL24 (2.1.2)**

587
588 Laboratory name: _____ Test name: _____ Test no. : _____
589

590 1. LabCyte CORNEA-MODEL24
591

592 Date received: _____
593

594 Lot no.: _____
595

596 Expiration date: _____
597

598
599 Accessories: Assay medium, 30mL (Lot no.: _____ Expiration date: _____)
600 24 well assay plate
601
602
603

<u>Note</u>

604

605

606 2. Assay medium
607
608

609 Date received: _____
610

611 Lot no.: _____
612

613 Expiration date: _____
614
615

<u>Note</u>

616

617

618 Date: _____ Operator: _____ Check date: _____ Study director: _____
619
620

621 Secretariat Check date: _____ Name: _____
622
623

624

625 **MDS 1-2:**
626 TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT
627 STEP-1 (3.2.1.1)

628 Laboratory name: _____ Test name: _____ Test no. : _____
629

630 1. Add distilled water (0.5 mL) to the wells of the 24-well assay plate.
631

632 To add distilled water (0.5 mL) Execution date/time: _____
633

634 2. Apply test chemicals to the wells of the 24-well assay plate.
635

636 3. Culture the 24-well assay plate in CO₂ incubator for 4 hours.
637

638 Time of WST-8 reaction started: _____
639

640 Time of WST-8 reaction completed: _____
641

642 4. Check the color of water.
643

644 5. Test chemical information and check list of coloring potential.
645
646
647

Test chemical	Physical state	Amount	Coloring	Test chemical.	Physical state	Amount	Coloring
PBS (NC)	LIQUID			Non treatment (NC)			
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	

648
649

Note

650 Date: _____ Operator: _____ Check date: _____ Study director: _____
651

652 Secretariat Check date: _____ Name: _____
653
654
655
656

657

658
659
660
661
662
663
664
665
666
667
668
669
670
671
672
673
674
675
676
677
678
679
680
681
682
683
684
685
686
687
688
689

MDS 1-3:
TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT
STEP-3 (3.2.2.1)

Laboratory name: _____ Test name: _____ Test no. : _____

- Preparation of WST-8 dilution medium

Warm EBSS for 30 minutes. Time/date: _____

CCK-8: (Lot no.: _____ Expiration date: _____)
EBSS: (Lot no.: _____ Expiration date: _____)
Volume _____ mL Time/date completed: _____
- Add WST-8 dilution medium (0.3mL) to the wells of the 24-well assay plate.

To add WST-8 dilution medium (0.3mL) Time/date executed: _____
- Apply test chemicals to the wells of the 24-well assay plate.
- Culture the 24-well assay plate in CO₂ incubator for 4 hours.

Time of WST-8 reaction started: _____
Time of WST-8 reaction completed: _____
- Check the color of WST-8 medium.
- Test chemical information and checked list of WST-8 assay interfere.


Test chemical	Physical state	Amount	WST-8 assay interfere	Test chemical.	Physical state	Amount	WST-8 assay interfere
PBS (NC)	LIQUID			Non treatment (NC)			
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	
	LIQUID	50 µL			SOLID	mg	

Note

Date: _____ Operator: _____ Check date: _____ Study director: _____
Secretariat Check date: _____ Name: _____

690
691

692
693
694
695
696
697
698
699

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 33 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

700 **MDS 1-4:**
701 **PREPARATION OF FREEZE KILLED TISSUE (3.2.2)**

702
703 Laboratory name: _____ Test name: _____ Test no. : _____

704 1. Transfer LabCyte CORNEA-MODEL24 tissues to 50 mL tube or appropriate sterile container.

705
706
707
708 2. Freeze tissues in the -80°C deep-freezer for 30 minutes (1st freezing).

709 Store for 30 minutes. Time/date: _____

710
711
712
713 3. Thaw tissues in the 37°C incubator for 15 minutes.

714 Store for 15 minutes. Time/date: _____

715
716
717
718 4. Freeze tissues in the -80°C deep-freezer for more than 30 minutes (2nd freezing).

719 Store for more 30minutes. Time/date: _____

720
721
722
723 5. Just before using, thaw tissues in the 37°C incubator for 15 minutes.

724 Store for 15 minutes. Time/date: _____

725
726
727
728
729
730
731
732
733

Note

734
735 Date: _____ Operator: _____ Check date: _____ Study director: _____

736
737
738 Secretariat Check date: _____ Name: _____

739
740
741

780 **MDS 3-1(LIQUID):**
781 **APPLICATION OF LIQUID TEST CHEMICALS, RINSING AND POST-INCUBATION (3.3.2)**

782 Laboratory name: _____ Test name: _____ Test no. : _____

783
784 1. Warm up the assay medium and add 0.5mL of the assay medium to the wells of the 3rd row on the
785 24-well assay plate for LIQUID.

786
787 Assay medium: (Lot no.: _____ Expiration date: _____)

788
789 Warm for 30 minutes. Time/date: _____
790 Add 0.5mL of assay medium to each well Time/date: _____ Number of plate: _____

791
792 2. Apply test chemicals to the LabCyte CORNEA-MODEL24.

793
794 Time/date execution started: _____

795
796 3. LIQUID test chemical information

Test chemical.	Lot no.	Physical state	Test chemical vol.	Time of application	Time of rinsing	Exposure period (1 minute)
PBS (Negative control)		LIQUID	50 µL	1	:	<input type="checkbox"/>
				2	:	<input type="checkbox"/>
				3	:	<input type="checkbox"/>
Ethanol (Positive control)		LIQUID	50 µL	1	:	<input type="checkbox"/>
				2	:	<input type="checkbox"/>
				3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	1	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	2	:	<input type="checkbox"/>
		LIQUID, viscous	50 µL	3	:	<input type="checkbox"/>

797
798 4. After exposure to test chemical for 1 minute, wash out the LabCyte CORNEA-MODEL24 and
799 transfer the culture inserts to the 3rd row on the 24-well assay plate.

800
801 PBS: (Lot no.: _____ Expiration date: _____)


802 Hit PBS stream on the tissue surface directly.

803
804 Confirm that there are no bubbles under the cell culture insert. Time/date completed: _____

805
806 5. Culture LabCyte CORNEA-MODEL24 in CO₂ incubator for 24 hours.

807
808 Time/date post-incubation started: _____

809 Time/date post-incubation completed: _____

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 36 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

810

<u>Note</u>

811
812
813
814
815
816
817
818
819
820
821
822
823

Date: _____ Operator: _____ Check date: _____ Study director: _____

Secretariat Check date: _____ Name: _____

MDS 3-2(SOLID):

APPLICATION OF SOLID TEST CHEMICALS AND RINSING (3.3.3)

Laboratory name: _____ Test name: _____ Test no. _____

1. Warm up the assay medium and add 0.5mL of the assay medium to the wells of 2nd and 3rd row on the 24-well assay plate for SOLID.

Assay medium: (Lot no.: _____ Expiration date: _____)

Warm for 30 minutes. Time/date: _____

Add 0.5mL of assay medium to each well Time/date: _____ Number of plate: _____

2. Apply test chemicals to the LabCyte CORNEA-MODEL24 and transfer the culture inserts to the 2nd row on the 24-well assay plate.

Time/date started: _____ Time/date completed: _____

3. SOLID test chemical information

Test chemical	Lot no.	Physical state	Crush and grind	Test chemical amount.	Time of application	Exposure period (24hours)
Non treatment (Negative control)					:	<input type="checkbox"/>
Lauric acid (Positive control)		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>
		SOLID	<input type="checkbox"/>	(mg, mg, mg)	:	<input type="checkbox"/>

4. Culture LabCyte CORNEA-MODEL24 in CO₂ incubator for 24 hours.

Time/date exposure started: _____


Time/date exposure completed: _____

5. After exposure to test chemical, wash out the LabCyte CORNEA-MODEL24 and transfer the culture inserts to the 3rd row on the 24-well assay plate.

PBS: (Lot no.: _____ Expiration date: _____)

Hit PBS stream on the tissue surface directly.

Time/date started: _____ Time/date completed: _____


Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 38 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

Note

860
861
862
863
864
865

Date: _____ Operator: _____ Check date: _____ Study director: _____

Secretariat Check date: _____ Name: _____

Version 2.5.6	<i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24	Page 39 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.

866 **MDS 4-1:**
867 **WST-8 ASSAY (3.3.4.1, 3.3.4.2)**

868
869 Laboratory name: _____ Test name: _____ Test no. : _____
870

871 1. Preparation of WST-8 dilution medium
872

873 Warm EBSS for 30 minutes. Time/date: _____
874
875 CCK-8: (Lot no.: _____ Expiration date: _____)
876 EBSS: (Lot no.: _____ Expiration date: _____)
877 Volume _____ mL Time/date completed: _____
878

879 2. Add WST-8 dilution medium (0.3mL) to the wells in the 4th row on the 24-well assay plate.
880
881 To add WST-8 dilution medium (0.3mL) Time/date executed: _____
882

883
884 3. After the operation, the blotted tissue transfer to wells of 4th row of 24-well assay plate.
885
886 Time/date started: _____ Time/date completed: _____
887
888 Confirm that there are no bubbles under the cell culture insert.

889
890 4. Culture LabCyte CORNEA-MODEL24 in CO₂ incubator for 4 hours.
891
892
893 Time of WST-8 reaction started: _____
894
895 Time of WST-8 reaction completed: _____
896
897
898
899
900
901
902
903
904

<u>Note</u>

905
906
907
908
909 Date: _____ Operator: _____ Check date: _____ Study director: _____
910
911
912 Secretariat Check date: _____ Name: _____
913
914
915

916 **MDS 4-2:**
917 SAMPLING OF REACTING WST-8 SOLUTION (3.3.4.3) AND MEASUREMENT (3.3.4.4)
918

919 Laboratory name: _____ Test name: _____ Test no. : _____

920 1. Reacting WST-8 solution (200 µL) is transferred to each well on the 96-well plate.

921 Transfer to the 96-well plate.

922 Time/date executed: _____

923 Sample location on 96-well plate.

	<u>LIQUID</u>						<u>SOLID</u>					
	1	2	3	4	5	6	7	8	9	10	11	12
A	blank											
B	PBS-1	PBS-2	PBS-3	Ethanol-1	Ethanol-2	Ethanol-3	Non-treatment-1	Non-treatment-2	Non-treatment-3	Lauric acid-1	Lauric acid-2	Lauric acid-3
C												
D												
E												
F												
G												
H												

929 2. Analyze extract OD at 450nm and 650nm, and calculate the OD(450nm-650nm).
930
931

932 Analyze OD at 450nm and 650nm.

933 Calculate the OD (450nm-650nm).

934 Calculate cell viability and SD.

935 Cell viability and SD are recorded on a separate data sheet.

936 The data sheet is attached to the back of this sheet.


937 Check for input errors.

938 Time/date executed: _____

Note

943 Date: _____ Operator: _____ Check date: _____ Study director: _____


944 Secretariat Check date: _____ Name: _____

Version 2.5.6	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 41 of 44
February, 2017		 Japan Tissue Engineering Co., Ltd.


949 **REVISION HISTORY**

950


Rev.	CONTENT	Date Revised
Ver.1.1	1) First version	Aug., 2012
Ver.1.2	1) Revised clerical error.	Aug., 2012
Ver.1.3	1) Revised clerical error.	Sep., 2012
Ver.1.3PS	<p>Added the following supplementary explanation.</p> 1) Added more detail explanation about the conditions of WST-8 reaction. 2) Added the formula of SD.	Sep., 2012
Ver.2.1	1) Revised clerical error. 2) In the section 1.2 “BACKGROUND”, changed statement about animal testing (Draize test). 3) In the section 3.3. “TEST METHOD”, explained the washing protocol of the <u>LIQUID</u> and <u>SOLID</u> test chemical more briefly. 4) In the section 3.3. “TEST METHOD”, changed the WST-8 dilution rate with PBS from 1:10 to 1:5. 5) In the section 3.3. “TEST METHOD”, changed the reaction period of WST-8 from 5 hours to 4 hours. 6) In the section 3.3. “TEST METHOD”, changed the condition of WST-8 reaction from shaking to standing. 7) In the section 3.3. “TEST METHOD”, changed the application-amount of <u>SOLID</u> chemicals from 50mg to 10mg. 8) As assay acceptance criteria, added that SD (test chemicals) of tissue viability of 3 identically treated replicates $\leq 20\%$. 9) In the prediction model of this EIT, changed the cut-off value of the mean viability from 50% to 40%.	May, 2013

Version 2.5.6	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 42 of 44
February, 2017		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>

Ver.2.2	<ol style="list-style-type: none"> 1) Revised clerical error. 2) In the section 3.3. "TEST METHOD", additionally explained that the solid test chemical is crush and grind in a mortar with pestle if necessary. 3) In the section 3.3.2 "APPLICATION OF LIQUID TEST CHEMICALS AND RINSING (DAY 0~1)", explained temperature condition (room temperature) at chemical application. 	Sep., 2013
Ver.2.3	<ol style="list-style-type: none"> 1) Revised clerical error. 2) In the section 3.3. "TEST METHOD", explained about a blank preparation at WST-8 assay. 3) In the section 3.3. "TEST METHOD", changed the condition of WST-8 reaction from standing to shaking. 5) In the section 3.3. "TEST METHOD", explained about a blank preparation at WST-8 assay. 	Feb., 2014
Ver.2.3.1	<ol style="list-style-type: none"> 1) Revised clerical error. 2) About the prediction result of eye irritation, changed the classification from irritation/no irritation to GHS classification. 	Mar., 2014
Ver.2.3.2MTT	<ol style="list-style-type: none"> 1) At the analysis of cell viability, changed the assay method from WST-8 assay to MTT assay. 	Jul., 2014
Ver.2.4.1	<ol style="list-style-type: none"> 1) At the analysis of cell viability, changed the assay method from MTT assay to WST-8 assay. 2) In the section 3.3. "TEST METHOD", changed the condition of WST-8 reaction from shaking to standing. 3) In the section 3.3. "TEST METHOD", changed the dilution solution of WST-8 reaction from PBS to EBSS. 4) In the section 3.3. "TEST METHOD", changed the reaction period of WST-8 from 5 hours to 4 hours. 5) The standard of the additional testing was mentioned about a 	Jan., 2015

Version 2.5.6	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 43 of 44
February, 2017		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>

	<p>borderline result.</p> <p>9) In the prediction model of this EIT, changed the cut-off value of the mean viability from 40% to 50%.</p>	
Ver.2.4.2	1) Changed the assay acceptance criterion of negative control from 0.5 < and <2.0 to 0.6 < and <1.5.	Mar., 2015
Ver.2.4.2m	1) The judgement of the chemical which is a result of the borderline was added.	Jul., 2015
Ver.2.5.1	<p>1) Revised clerical error.</p> <p>2) In the prediction model of this EIT, changed the cut-off value of the mean viability from 50% to 40%.</p> <p>3) Changed the judgement of equivocal results.</p> <p>4) Changed the assay acceptance criterion of negative control from 0.6 < and <1.5 to 0.5 < and <1.3.</p> <p>5) In the section 3.2. "TEST FOR DETECTING CHEMICALS THAT INTERFERE WITH WST-8 ENDPOINT", revised description has clearly.</p> <p>6) Added MDS1-2 and MDS 1-3.</p>	Aug., 2015
Ver.2.5.1m	<p>1) Revised clerical error.</p> <p>2) Changed allocation for a 96-well plate of pattern which consists of freeze killed tissue (Fig.12B).</p>	Sep., 2015
Ver. 2.5.1mr	1) Revised Fig.2B	Sep.,2015
Ver. 2.5.2	<p>1) In the section 3.3. "TEST METHOD", explained the washing protocol of the <u>LIQUID</u> and <u>SOLID</u> test chemical more briefly.</p> <p>2) In the MDS 3-1 and the MDS 3-2, added the check box about the washing procedure.</p> <p>3) In the MDS 3-2, added the check box about crush and grind of</p>	Sep.,2015

Version 2.5.6	<p style="text-align: center;"><i>IN VITRO</i> EYE IRRITATION TEST USING HUMAN CORNEAL TISSUE MODEL: LabCyte CORNEA-MODEL24</p>	Page 44 of 44
February, 2017		 <p style="text-align: center;">Japan Tissue Engineering Co., Ltd.</p>

	test chemicals.	
Ver. 2.5.4	<ol style="list-style-type: none"> 1) In the section 3.3. "TEST METHOD", added attention point of the washing procedure. 2) In the section 4.3 "ASSESSMENT FLOWCHART", Change numbers of test run from three independent run to single run. 	June, 2016
Ver. 2.5.5	<ol style="list-style-type: none"> 1) Revised clerical error. 2) As assay acceptance criteria, changed that SD (negative control, positive control and test chemicals) of tissue viability of 3 identically treated replicates from $\leq 20\%$ to $\leq 18\%$. 	October, 2016
Ver.2.5.6	<ol style="list-style-type: none"> 1) Revised clerical error. 2) Detection protocol of coloring interference is changed from using WST-8 medium to distilled water and correction of coloring interference is changed from using freeze-killed tissue to using living tissue without WST-8 reaction. 	February, 2017

951

952