

Appendix B1

Draft ICCVAM Recommended Protocol for Future Studies Using the Cryopreserved Whole Blood (Cryo WB)/Interleukin-1 β (IL-1 β) Test Method

[This Page Intentionally Left Blank]

1 **Draft ICCVAM Recommended Protocol for Future Studies Using the Cryopreserved**
2 **Whole Blood (Cryo WB)/Interleukin-1 β (IL-1 β) Test Method**

3

4

PREFACE

5

6 This proposed protocol for the detection of pyrogenicity is based on information obtained
7 from 1) The European Centre for the Validation of Alternative Methods (ECVAM) Cryo
8 WB/IL-1 β Background Review Document (BRD) presented in Appendix A of the draft
9 Interagency Coordination Committee on the Validation of Alternative Methods (ICCVAM)
10 Pyrogenicity Test Method BRD, which includes ECVAM Standard Operating Procedures
11 (SOPs) for the Cryo WB/IL-1 β test method, and 2) Information provided to the National
12 Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative
13 Toxicological Methods (NICEATM) by Dr. Thomas Hartung, Head of ECVAM. The
14 ECVAM SOPs are based on the Cryo WB/IL-1 β methodology first described by Schindler et
15 al. (2004). A table of comparison between the draft ICCVAM recommended protocol and the
16 ECVAM SOP is provided in **Table 1**. Future studies using the Cryo WB/IL-1 β test method
17 may include further characterization of the usefulness or limitations of the assay for
18 regulatory decision-making. Users should be aware that the proposed test method protocol
19 might be revised based on additional optimization and/or validation studies. ICCVAM
20 recommends that test method users routinely consult the ICCVAM/NICEATM website
21 (<http://iccvam.niehs.nih.gov/>) to ensure that the most current test method protocol is used.

22

22 **Table 1 Comparison of Draft ICCVAM Recommended Test Method Protocol**
 23 **with the ECVAM SOP for the Cryopreserved Whole Blood (Cryo**
 24 **WB)/Interleukin-1 β (IL-1 β) Pyrogen Test**

Protocol Component	Draft ICCVAM Protocol	ECVAM SOP
Test Substance	Test neat or at minimal dilution that produces 50% to 200% of 1 EU/mL EC	Test at MVD
Incubation Plate (number of control or test groups at n=4 each)	NSC (1)	Same as ICCVAM protocol
	EC (5)	EC (2)
	TS (14)	TS (3) x EC (5) spikes
	PPC ¹ (0)	PPC (3)
	NPC ¹ (0)	NPC (3)
ELISA Plate	Includes seven point IL-1 β SC and blank in duplicate	Same as ICCVAM protocol
Decision Criteria for Interference	0.5 x Median OD ₄₅₀ ² of 1 EU/mL EC <2x Median OD ₄₅₀ of 1 EU/mL EC	Mean OD ₄₅₀ ² of PPC \geq 1.6X Mean OD ₄₅₀ of NPC
Assay Acceptability Criteria	Mean OD ₄₅₀ of 0.5 EU/mL EC \geq 1.6X Mean OD ₄₅₀ of NSC	Same as ICCVAM protocol
	Mean OD ₄₅₀ of PPC is 50% to 200% of 0.5 EU/mL EC	Same as ICCVAM protocol
	Mean OD ₄₅₀ of NSC \leq 0.15 ³	Mean OD ₄₅₀ of NSC \leq 0.10
	Not included	If one OD ₄₅₀ of the 1.0 EU/mL EC >Max, the ELISA may be repeated at reduced incubation time
	EC SC produces OD ₄₅₀ values that ascend in a sigmoidal concentration response	Same as ICCVAM protocol
	Quadratic function of IL-1 β SC r ² \geq 0.95 ³	Same as ICCVAM protocol
Decision Criteria for Pyrogenicity	OD ₄₅₀ Test sample \geq OD ₄₅₀ 0.5 EU/mL EC	Same as ICCVAM protocol

25 Abbreviations: EC = Endotoxin control; MVD = Maximum Valid Dilution; NPC = Negative Product Control;
 26 NSC = Normal saline control;

27 PPC = Positive Product Control; SC = Standard curve; TS = Test substance

28 ¹ PPC and NPC are evaluated during the interference test.

29 ² Median or mean OD₄₅₀ values are corrected (i.e., reference filter reading, if applicable, and NSC are
 30 subtracted).

31 ³ Criteria originated from PBMC SOP.
 32
 33
 34

34 **1.0 PURPOSE AND APPLICABILITY**

35 The purpose of this protocol is to describe the procedures used to evaluate the presence of a
36 pyrogen (i.e., Gram-negative endotoxin) in a test substance. The presence of Gram-negative
37 endotoxin is detected by its ability to induce cytokine IL-1 β release from monocytoid cells in
38 human Cryo WB. The quantity of IL-1 β released is obtained using an enzyme-linked
39 immunosorbent assay (ELISA) that includes monoclonal or polyclonal antibodies specific for
40 IL-1 β . Release of IL-1 β is measured by incubation of Cryo WB with test substances or
41 controls (i.e., positive, negative). The amount of pyrogen present is determined by comparing
42 the values of endotoxin equivalents produced by cells exposed to the test substance to those
43 exposed to an internationally-harmonized Reference Standard Endotoxin (RSE)¹ or an
44 equivalent standard expressed in Endotoxin Units (EU)/mL. Based on a rabbit threshold
45 pyrogen dose of 0.5 EU/mL, which was established in a retrospective evaluation of rabbit
46 pyrogen test (RPT) data, a test substance is considered pyrogenic if it induces a level of IL-
47 1 β release equal to or greater than that induced by 0.5 EU/mL of endotoxin.

48 The focus of this protocol is on the use of the Cryo WB/IL-1 β test method, specifically for
49 the detection of Gram-negative endotoxin in parenteral pharmaceuticals. The relevance and
50 reliability for non-endotoxin pyrogens (e.g., lipoteichoic acid) has not been demonstrated in a
51 formal validation study.

52 **2.0 SAFETY AND OPERATING PRECAUTIONS**

53 All procedures for procurement of eligible blood donors and blood donations should follow
54 the regulations and procedures set forth by institutional guidelines for utilization of human
55 substances, which include but are not limited to blood, tissues, and tissue fluids. Standard
56 laboratory precautions are recommended including the use of laboratory coats, eye
57 protection, and gloves. If necessary, additional precautions required for specific study
58 substances or hazardous chemicals will be identified in the Material Safety Data Sheet
59 (MSDS).

¹ RSEs are internationally-harmonized reference standards (e.g., WHO *E. coli* Lipopolysaccharide [LPS] 94/580 [0113:H10:K-]; USP RSE Lot G3E069; FDA Lot EC-6). Equivalent endotoxins include commercially available *E. coli*-derived LPS Control Standard Endotoxin (CSE) or other *E. coli* LPS preparations that have been calibrated with an appropriate RSE.

60 The stop solution used in the ELISA kit is acidic and corrosive and should be handled with
61 the proper personal protective devices. If this reagent comes into contact with skin or eyes,
62 wash thoroughly with water. Seek medical attention, if necessary.

63 Tetramethylbenzidine (TMB) solution contains a hydrogen peroxide substrate and 3, 3', 5, 5'-
64 TMB. This reagent is a strong oxidizing agent and a suspected mutagen. Appropriate
65 personal protection should be used to prevent bodily contact.

66 Bacterial endotoxin is a toxic agent (i.e., can induce sepsis, shock, vascular damage,
67 antigenic response) and should be handled with care. Skin cuts should be covered and
68 appropriate personal protective devices should be worn. In case of contact with endotoxin,
69 immediately flush eyes or skin with water for at least 15 min. If inhaled, remove the affected
70 individual from the area and provide oxygen and/or artificial respiration as needed. Skin
71 absorption, ingestion, or inhalation may produce fever, headache, and hypotension.

72 **3.0 MATERIALS, EQUIPMENT, AND SUPPLIES**

73 **3.1 Source of Cells**

74 Leukocytes from WB are the source of cells for cytokine production in the Cryo WB/IL-1 β
75 test method (Hartung and Wendel, 1999; Schindler et al., 2004, 2006). WB is obtained from
76 healthy human volunteers who have provided their consent according to established
77 institutional guidelines. Volunteers are expected not to have taken any drugs (e.g.,
78 prescription drugs, recreational drugs, herbal drugs) and to have been free from illness for at
79 least two weeks prior to donation.

80 The WB is processed and cryopreserved using either the Konstanz method developed at the
81 University of Konstanz (Schindler et al., 2004) or the PEI method developed at the Paul
82 Ehrlich Institute (Schindler et al., 2006).

83 **3.2 Equipment and Supplies**

84 For all steps in the protocol, excluding the ELISA procedure, the materials that will be in
85 close contact with samples and/or blood cells (e.g., pipet tips, containers, solutions) should be
86 sterile and free from detectable pyrogens.

- 87 3.2.1 Blood Incubation
- 88 3.2.1.1 *Equipment*
- 89 • Centrifuge
- 90 • Hood; Bio-safety, laminar flow (recommended)
- 91 • Incubator; cell culture (37 \pm 1 $^{\circ}$ C + 5% CO₂)
- 92 • Pipetter; multichannel (8- or 12-channel)
- 93 • Pipettors; single-channel adjustable (20 and 200 μ L)
- 94 • Repeating pipetter
- 95 • Vortex mixer
- 96 3.2.1.2 *Consumables*
- 97 • Centrifuge tubes; nonpyrogenic, polystyrene (15 and 50 mL)
- 98 • Combitips; repeating pipetter (1.0 and 2.5 mL)
- 99 • Needle set; Sarstedt multify, pyrogen-free, 19 mm, 21 gauge for S-Monovette
- 100 • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 101 • Reaction tubes; polystyrene (1.5 mL)
- 102 • Reservoirs; fluid
- 103 • RPMI-1640 cell culture medium
- 104 • Tips; pipetter, sterile, pyrogen-free (20 and 200 μ L)
- 105 • Tubes; Sarstedt S-Monovette, 7.5 mL, heparinized for blood collection
- 106 3.2.2 ELISA
- 107 3.2.2.1 *Equipment*
- 108 • Microplate mixer
- 109 • Microplate reader (450 nm with an optional reference filter in the range of
- 110 600-690 nm)
- 111 • Microplate washer (optional)

- 112
- Multichannel pipetter

113 3.2.2.2 *Consumables*

- 114
- Container; storage, plastic
- 115
- Deionized water; nonsterile
- 116
- Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 117
- Pyrogen-free water (PFW)
- 118
- Reservoirs; fluid
- 119
- Tips; pipetter, nonsterile
- 120
- Tubes; polystyrene (12mL)

121 3.2.2.3 *ELISA Kit*

122 An ELISA that measures IL-1 β release from Cryo WB is used. A variety of IL-1 β ELISA
123 kits are commercially available and the IL-1 β ELISA procedure outlined in this protocol is
124 intended to serve as an example for using an ELISA kit. If the user prefers to prepare an in-
125 house ELISA, then additional reagents would be required. The IL-1 β ELISA should be
126 calibrated using an international reference standard (e.g., WHO 86/680) prior to use. The IL-
127 1 β cytokine assay kits do not provide the RSE or endotoxin equivalent; therefore, this
128 reagent must be purchased separately. Results obtained using these products are subject to
129 the assay acceptability and decision criteria described in **Sections 8.0** and **9.0**. IL-1 β ELISA
130 kit components may include the following:

- 131
- ELISA plates coated with anti-human IL-1 β capture antibody; monoclonal or
- 132
- polyclonal
- 133
- Buffered wash solution
- 134
- Dilution buffer
- 135
- Enzyme-labeled detection antibody
- 136
- Human IL-1 β reference standard
- 137
- Pyrogen-free saline (PFS)

- 138 • Stop solution
- 139 • TMB/substrate solution

140 **3.3 Chemicals**

- 141 • Endotoxin (e.g., WHO *E. coli* LPS 2nd International Standard 94/580; USP
- 142 RSE *E. coli* LPS Lot G3E069; USP RSE *E. coli* Lot G; FDA *E. coli* Lot EC6)

143 **3.4 Solutions**

- 144 • RPMI-1640 cell culture medium

145 **4.0 ASSAY PREPARATION**

146 All test substances, endotoxin, and endotoxin-spiked solutions should be stored at 4°C.

147 **4.1 Endotoxin Standard Curve**

148 An internationally harmonized RSE or equivalent is used to generate the endotoxin standard
149 curve. The use of any other *E. coli* LPS requires calibration against a RSE using the Cryo
150 WB/IL-1 β test method.

151 A standard endotoxin curve consisting of a Normal Saline Control (NSC) and five RSE
152 concentrations (0.25, 0.50, 1.0, 2.5, and 5.0 EU/mL) are included in the incubation step (refer
153 to **Table 4-1**) and then transferred to the ELISA plate. To prepare the endotoxin standard
154 curve, first obtain a 2000 EU/mL stock solution by addition of PFW to the lyophilized
155 content of the stock vial by following the instructions provided by the manufacturer (e.g., for
156 a vial containing 10,000 EU, 5 mL of PFW is added). To reconstitute the endotoxin, the stock
157 vial should be vortexed vigorously for at least 30 min or sonicated in a bath sonicator for 5
158 min. The stock solution is stable for 14 days when stored at 2 to 8°C. An endotoxin standard
159 curve is prepared by making serial dilutions of the stock solution in PFS as described in
160 **Table 4-1**.

161

161

162 **Table 4-1 Preparation of Endotoxin Standard Curve**

Stock Endotoxin EU/mL	μ L of Stock Endotoxin	μ L of PFS	Endotoxin Concentration EU/mL
2000 ^{1,2}	50	1950	50 ³
50	100	900	5.0
5.0	500	500	2.5
2.5	400	600	1.0
1.0	500	500	0.50
0.50	500	500	0.25
0	0	1000	0

163 Abbreviations: EU = Endotoxin units; PFS = Pyrogen-free saline

164 Each stock tube should be vortexed vigorously prior to its use to make the subsequent dilution.

165 ¹ A 2000 EU/mL stock solution of endotoxin is prepared according to the manufacturer's instructions.166 ² The stock solution of USP RSE may be stored in aliquots and kept at -20°C for up to 6 months. Do not store
167 the endotoxin at -80°C.168 ³ This concentration is not used in the assay.

169

170 **4.2 Test Substances**

171 Liquid test substances should be tested neat or, if interference is detected (see **Section 4.2.1**),
172 diluted in PFS. Solid test substances should be prepared as solutions in PFS or, if insoluble in
173 saline, dissolved in dimethylsulfoxide (DMSO) then diluted up to 0.5% (v/v) with PFS,
174 provided that this concentration does not interfere with the assay. The test substances should
175 be vortexed vigorously for at least 30 min or sonicated in a bath sonicator for 5 min.

176 **4.2.1 Interference Testing**

177 Interference testing must be carried out on any test sample for which no interference
178 information is available. The purpose of the interference test is to determine the lowest
179 dilution (i.e., highest concentration) of a test substance from which an endotoxin spike can be
180 detected (i.e., based on the decision criteria described in **4.2.1.2**). However, to ensure a valid
181 test, a test substance should not be diluted beyond its Maximum Valid Dilution (MVD).

182 For many marketed products, values for the MVD and the Endotoxin Limit Concentration
183 (ELC) are published in the U.S. Pharmacopeia, the European Pharmacopoeia, and/or Food
184 and Drug Administration (FDA) guidelines. However if one or both of these values are not
185 available, then calculation of the MVD is dependent on the ELC (see **Section 12.3**). If
186 unknown, the ELC can be approximated by dividing the maximum hourly dose of the

187 product by the hourly dose received per patient. For example, if a product is used at an
188 hourly dose of 100 mg per patient, then the ELC would be 350 EU/100 mg or 3.5 EU/mg.

189 4.2.1.1 Reference Endotoxin for Spiking Test Substances

190 The WHO-LPS 94/580 [*E. coli* O113:H10:K-] or equivalent international RSE (e.g., USP
191 G3E069, FDA EC-6) is recommended for preparation of the endotoxin control (EC). If a
192 different *E. coli* LPS is used and the potency relative to the RSE is not provided, then each
193 lot must be calibrated against the RSE in the Cryo WB/IL-1 β test method. For interference
194 testing, an endotoxin standard curve (see **Section 4.1**) should be included on each plate.

195 4.2.1.2 Spiking Test Substances with Endotoxin

196 For interference testing, non-spiked and endotoxin-spiked test substances are prepared in
197 microplate wells (n=4 replicates) and an *in vitro* pyrogen test is performed. Either RPMI or a
198 fixed concentration (a concentration selected from the middle of the EC standard curve) of
199 the RSE (i.e., 1 EU/mL) in RPMI is added to the test substance in serial two-fold dilutions.
200 An illustrative example of endotoxin spiking solutions is shown in **Table 4-2**. For non-spiked
201 solutions, 200 μ L of RPMI is added to a well followed by 20 μ L of the test substance (neat or
202 at serial dilution) and 20 μ L of WB. Endotoxin-spiked solutions are prepared by adding 180
203 μ L of RPMI to each well followed by 20 μ L of the test substance (neat or at serial dilution)
204 and 20 μ L of WB. Then, 20 μ L of a 1 EU/mL solution of endotoxin in RPMI is added and
205 the well contents are mixed (see example presented in **Table 4-2**).

206

206

207 **Table 4-2 Preparation of Endotoxin-Spiked and Non-spiked Solutions for**
 208 **Determination of Test Substance Interference in the Incubation and**
 209 **ELISA Test Systems**

Sample Addition	Spiked	Non-spiked
	$\mu\text{L}/\text{well}^1$	
RPMI	180	200
Endotoxin spike solution ²	20	0
Test substance (neat and each serial dilution)	20	20
Cryo WB	20	20
Total ³	240	240

210 Abbreviations: Cryo WB = Cryopreserved whole blood

211 ¹ n=4 replicates each

212 ² Endotoxin concentration is 1.0 EU/mL in RPMI.

213 ³ A total volume of 240 μL per well is used for the incubation.

214

215

216 The lowest dilution of the test substance that yields an endotoxin spike recovery of 50% to
 217 200% in the pyrogen test is determined. The optical density (OD) values of the endotoxin-
 218 spiked and non-spiked test substances are calibrated against the endotoxin calibration curve.
 219 The resulting EU value of the non-spiked test substance is subtracted from the corresponding
 220 EU value of the endotoxin-spiked test substance at each dilution. The % recovery for each
 221 sample dilution is then determined from the endotoxin spike solution concentration set to
 222 100%. For example, consider the following interference test results in **Table 4-3**:

223 **Table 4-3 Example of Interference Data Used to Determine Sample Dilution for**
 224 **Assay**

Sample Dilution	% Recovery of Endotoxin Control
None	25
1:2	49
1:4	90
1:8	110

225

226

227 Based on these results, the dilution of the test substance used in the *in vitro* pyrogen test
 228 would be 1:4 (i.e., the lowest dilution between 50% and 200% of the 1 EU/mL EC).

229 4.2.2 Interference with ELISA System

230 If the data obtained from the experiment in **Section 4.2.1** suggests the presence of
231 interference, then a subsequent experiment similar to that described in **Section 4.2.1** would
232 need to be performed to confirm that the test substance(s) does not directly interfere with the
233 ELISA. For this experiment, an ELISA would be performed in the absence of Cryo WB.

234 **5.0 CONTROLS**

235 **5.1 Negative Control**

236 A negative control (e.g., PFS is added instead of the test sample) is included in each
237 experiment in order to detect nonspecific changes in the test system, as well as to provide a
238 baseline for the assay endpoints.

239 **5.2 Solvent Control**

240 Solvent controls are recommended to demonstrate that the solvent is not interfering with the
241 test system when solvents other than PFS are used to dissolve test substances.

242 **5.3 Positive Control**

243 An internationally standardized EC (e.g., WHO 94/580; USP G3E069; 0.5 EU/mL) is
244 included in each experiment to verify that an appropriate response is induced.

245 **5.4 Benchmark Controls**

246 Benchmark controls may be used to demonstrate that the test method is functioning properly,
247 or to evaluate the relative pyrogenic potential of chemicals (e.g., parenteral pharmaceuticals,
248 medical device eluates) of a specific class or a specific range of responses, or for evaluating
249 the relative pyrogenic potential of a test substance. Appropriate benchmark controls should
250 have the following properties:

- 251 • consistent and reliable source(s) for the chemicals (e.g., parenteral
252 pharmaceuticals, medical device eluates)
- 253 • structural and functional similarities to the class of substance being tested
- 254 • known physical/chemical characteristics

- 255 • supporting data on known effects in animal models
- 256 • known potency in the range of response

257 **5.5 Positive Product Control (PPC)**

258 The PPC is a test substance diluted to a level that does not interfere with the test method and
259 does not exceed the MVD. The PPC is obtained by spiking a test substance with a known EC
260 (e.g., 1 EU/mL) and demonstrating that 50% to 200% of the EC is recovered.

261 **5.6 Negative Product Control (NPC)**

262 The NPC is the test substance diluted to the MVD and then spiked with PFS. It is the
263 negative control for the PPC.

264 **6.0 EXPERIMENTAL DESIGN**

265 **6.1 Incubation with Test Samples and Measurement of IL-1 β Release**

266 6.1.1 Collection of Human Blood

267 WB is obtained from healthy human volunteers who have provided their consent according
268 to established institutional guidelines. Volunteers are expected not to have taken any drugs
269 and to have been free from illness for at least two weeks prior to donation. The criteria for
270 rejection of data from donors that are low responders or that are suspect due to veracity of
271 health information is addressed in **Section 8.0**.

272 Fresh WB is drawn by venipuncture using a multily needle set and collected in heparinized
273 tubes (15 IU/mL lithium heparin). All components of the blood collection system (i.e.,
274 syringes, tubes, connecting lines) must be sterile and pyrogen-free.

275 6.1.2 Cryopreservation Procedure

276 Two methods are available for cryopreservation of blood 1) The PEI method developed at the
277 Paul Ehrlich Institute and 2) The Konstanz method developed at the University of Konstanz.

278 6.1.3 PEI Method of Cryopreservation

279 In the PEI method (Schindler et al., 2006), heparinized WB pooled from five donors is frozen
280 at -80°C in a cryoprotective phosphate buffer (Soerensen's) containing 10% (v/v) pyrogen-
281 free, clinical-grade DMSO in cryotubes.

282 6.1.3.1 *Konstanz Method of Cryopreservation*

283 In the Konstanz method (Schindler et al., 2004), pyrogen-free, clinical grade DMSO is added
284 directly to the blood of individual donors at a final concentration of 10% (v/v). The blood is
285 then pooled and 1.2 mL aliquots are placed in cryotubes. The blood is frozen in a computer-
286 controlled freezer using several cycles of programmed freezing down to -120°C. Tubes of
287 blood are then removed from the instrument and placed in liquid nitrogen.

288 6.1.3.2 *Thawing Procedure*

289 Calculate the volume of Cryo WB needed to carry out the assay (20 μ L/well or 1.92 mL/96-
290 well plate) and remove a sufficient number of aliquots from the freezer. Place the tubes in an
291 incubator at 37 \pm 1°C and allow them to thaw for 15 min. In a laminar flow hood, unscrew the
292 caps and pool the Cryo WB in a centrifuge tube. Mix the tubes by gentle inversion. **Do Not**
293 **Vortex.**

294 6.1.4 Incubation Plate

295 Test substances are prepared at a level of dilution that did not show interference with the test
296 system or for which it is known that interference does not occur. Each incubation plate can
297 accommodate an endotoxin standard curve, a NSC, and 14 test samples (see **Table 6-1**).
298

298

299 **Table 6-1 Overview of Incubation Plate Preparation in the Cryo WB/IL-1 β Test**
 300 **Method (PEI Method)**

Number of Wells	Sample	RPMI	EC	Test Sample	Cryo WB ¹	Mix the samples; incubate overnight at 37 \pm 1 $^{\circ}$ C in a humidified atmosphere with 5% CO ₂ .	Mix the samples; immediately transfer to an ELISA plate ⁴ and run ELISA or store plate at -20 $^{\circ}$ C or -80 $^{\circ}$ C.
20 ²	EC	180	20	0	40		
4	NSC	220	0	0	20		
56 ³	Test samples (1-14)	200	0	20	20		

301 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; Cryo WB = Cryopreserved whole
 302 blood

303 ¹ For the Konstanz method of cryopreservation, 20 μ L of Cryo WB is used and the volume of RPMI is adjusted
 304 to 200 μ L.

305 ² Five EC concentrations (0.25, 0.50, 1.0, 2.5, 5.0 EU/mL) in quadruplicate.

306 ³ 14 test samples (n=4) per plate.

307 ⁴ An IL-1 β standard curve is prepared in Columns 11 and 12 on the ELISA plate. Therefore, 80 wells are
 308 available for test samples and controls on the incubation plate.

309

310 6.1.5 Incubation Assay for IL-1 β Release

311 Blood samples are prepared in a microtiter plate using a laminar flow hood. All consumables
 312 and solutions must be sterile and pyrogen-free. Each plate should be labeled appropriately
 313 with a permanent marker. An overview of the incubation plate preparation is shown in **Table**
 314 **6-1**. The incubation procedure is outlined below:

315 **Step 1.** Refer to the incubation plate template presented in **Table 6-2**.

316 **Step 2.** Using a pipetter, transfer either 200 or 180 μ L of RPMI into each well (for
 317 the Konstanz or PEI method of cryopreservation, respectively – refer to **Step 5**
 318 below).

319 **Step 3.** Transfer 20 μ L of test sample into the appropriate wells as indicated in the
 320 template.

321 **Step 4.** Transfer 20 μ L of the EC (standard curve) and the NSC controls in
 322 quadruplicate into the appropriate wells according to the template.

323 **Step 5.** Transfer either 20 or 40 μ L of Cryo WB (for the Konstanz or PEI method
324 of cryopreservation, respectively) into each well and mix by gently swirling the
325 plate.

326 **Step 6.** Mix the contents of the wells thoroughly by gently pipetting up and down
327 five times using a multichannel pipetter, changing the tips between each row in
328 order to avoid cross-contamination.

329 **Step 7.** Place the covered plate in a tissue culture incubator for 10 to 24 hr at
330 $37\pm 1^\circ\text{C}$ in a humidified atmosphere containing 5% CO_2 .

331 **Step 8.** If using the Konstanz method, freeze the plate at -20°C or -80°C until the
332 contents of the well are completely frozen and then, thaw the plate at RT or in a
333 water bath not exceeding 37°C .

334 **Step 9.** Prior to transferring the test samples onto the ELISA plate, mix the
335 contents of the wells by pipetting up and down three times using a multichannel
336 pipetter, changing the tips between each row in order to avoid cross-contamination.

337 *Note: The aliquots may be tested immediately in the ELISA or stored at -20°C or*
338 *-80°C for testing at a later time. After transfer to the ELISA plate, freeze the*
339 *remaining aliquots at -20°C or -80°C for subsequent experiments, if necessary (see*
340 *Assay Acceptability and Decision Criteria in **Sections 8.0 and 9.0**).*

341

341

342 **Table 6-2 Incubation Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
A	EC ¹ 5.0	EC 5.0	EC 5.0	EC 5.0	TS3	TS3	TS3	TS3	TS11	TS11	Void ³	Void
B	EC 2.5	EC 2.5	EC 2.5	EC 2.5	TS4	TS4	TS4	TS4	TS11	TS11	Void	Void
C	EC 1.0	EC 1.0	EC 1.0	EC 1.0	TS5	TS5	TS5	TS5	TS12	TS12	Void	Void
D	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS6	TS6	TS6	TS6	TS12	TS12	Void	Void
E	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS7	TS7	TS7	TS7	TS13	TS13	Void	Void
F	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	Void	Void
G	TS1 ²	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	Void	Void
H	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	Void	Void

343 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

344 ¹ EC value (e.g., EC 5.0) represents the endotoxin concentration in EU/mL.345 ² TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.346 ³ Columns 11 and 12 are reserved for the IL-1 β standard curve on the ELISA plate (see **Table 6-3**).

347

348

349 **6.2 ELISA to Measure IL-1 β Release**350 **6.2.1 IL-1 β Standard Curve**

351 An IL-1 β standard, supplied with the ELISA kit, is used. IL-1 β standards are typically
 352 supplied in lyophilized form and should be reconstituted according to the manufacturer's
 353 instructions. The stock solution should be diluted in RPMI to the following concentrations: 0,
 354 62.5, 125, 250, 500, 1000, 2000, and 4000 pg/mL. Each well on the ELISA plate will receive
 355 100 μ L of an IL-1 β blank or standard.

356 **6.2.2 ELISA**

357 The manufacturer's instructions provided with the ELISA kit should be followed and a
 358 typical experimental design is outlined below. If the user prefers to prepare an in-house
 359 ELISA, then appropriate modification and validation of these changes would be necessary.
 360 The ELISA should be carried out at room temperature (RT) and therefore all components
 361 must be at RT prior to use. Do *not* thaw frozen specimens by heating them in a water bath. A
 362 suggested ELISA plate template is shown in **Table 6-3**, which includes a five-point EC
 363 standard curve, a NSC, an eight-point IL-1 β standard curve (0 to 4000 pg/mL), and 14 test

364 substances in quadruplicate. The EC standard curve, the NSC, and the test sample
365 supernatants are transferred directly from the incubation plate. The IL-1 β standard curve is
366 prepared as described in **Section 6.2.1**. An overview of the ELISA plate preparation is shown
367 in **Table 6-4**.

368 **Step 1.** Add 100 μ L of enzyme-labeled detection antibody to each well.

369 **Step 2.** After pipetting up and down three times to mix the supernatant, transfer
370 100 μ L from each well of the Incubation Plate (A1-10; H1-10) to the ELISA plate.

371 **Step 3.** Add 100 μ L of each IL-1 β standard (0 to 4000 pg/mL) into the respective
372 wells on the ELISA plate.

373 **Step 4.** Cover the microtiter plate(s) with adhesive film and incubate for 90 min on
374 a microplate mixer at 350-400 rpm at 20 to 25°C.

375 **Step 5.** Decant and wash each well three times with 300 μ L Buffered Wash
376 Solution and then rinse three times with deionized water. Place the plates upside
377 down and tap to remove water.

378 **Step 6.** Add 200 μ L of TMB/Substrate Solution to each well and incubate at RT in
379 the dark for 15 min. If necessary, decrease the incubation time.

380 **Step 7.** Add 50 μ L of Stop Solution to each well.

381 **Step 8.** Tap the plate gently after the addition of Stop Solution to aid in mixing.

382 **Step 9.** Read the OD₄₅₀ within 15 min of adding the Stop Solution. Measurement
383 with a reference wavelength of 600-690 nm is recommended.

384

384

385 **Table 6-3 ELISA Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
A	EC ¹ 5.0	EC 5.0	EC 5.0	EC 5.0	TS3	TS3	TS3	TS3	TS11	TS11	IL-1 β ³ 0	IL-1 β 0
B	EC 2.5	EC 2.5	EC 2.5	EC 2.5	TS4	TS4	TS4	TS4	TS11	TS11	IL-1 β 62.5	IL-1 β 62.5
C	EC 1.0	EC 1.0	EC 1.0	EC 1.0	TS5	TS5	TS5	TS5	TS12	TS12	IL-1 β 125	IL-1 β 125
D	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS6	TS6	TS6	TS6	TS12	TS12	IL-1 β 250	IL-1 β 250
E	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS7	TS7	TS7	TS7	TS13	TS13	IL-1 β 500	IL-1 β 500
F	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	IL-1 β 1000	IL-1 β 1000
G	TS1 ²	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	IL-1 β 2000	IL-1 β 2000
H	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	IL-1 β 4000	IL-1 β 4000

386 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

387 ¹ EC value (e.g., EC 5.0) represents the endotoxin concentration in EU/mL.388 ² TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.389 ³ IL-1 β values in columns 11 and 12 are in pg/mL.

390

391

392 **Table 6-4 Overview of ELISA Procedure**

Enzyme-labeled Antibody (μL)	Material transfer from Incubation Plate (μL)	Incubate 90 min on a plate mixer at 350-400 rpm.	TMB/Substrate Solution (μL)	Incubate 15 min at RT in dark.	Stop Solution (μL)	Read optical density at 450 nm with a 600-690 nm wavelength reference filter.
100	100		200		50	

393 Abbreviations: RT = Room temperature

394

395 **7.0 EVALUATION OF TEST RESULTS**396 **7.1 OD Measurements**

397 The OD of each well is obtained by reading the samples in a standard microplate

398 spectrophotometer (i.e., plate reader) using a visible light wavelength of 450 nm (OD₄₅₀) with

399 a 600 to 690 nm reference filter (recommended). OD₄₅₀ values are used to determine assay
400 acceptability and in the decision criteria for the detection of endotoxin in a test substance (see
401 **Sections 8.0** and **9.0**).

402 **8.0 CRITERIA FOR AN ACCEPTABLE TEST**

403 Obtain the PPC and the corresponding NPC by interference testing of a test substance in the
404 presence and absence of a fixed quantity of endotoxin (i.e., 1.0 EU/mL) in quadruplicate. An
405 EC (five-point standard curve) and a NSC should be included in each experiment. An IL-6
406 standard curve should be included in each ELISA as shown in the template presented in
407 **Table 6-3**. An assay is considered acceptable only if the following minimum criteria are met:

- 408 • The quadratic function of the IL-1 β standard curve produces an $r \geq 0.95$ and
409 the OD₄₅₀ of the blank control is below 0.15.
- 410 • The endotoxin standard curve produces OD₄₅₀ values that ascend in a
411 sigmoidal concentration response.

412 **9.0 DATA INTERPRETATION/DECISION CRITERIA**

413 **9.1 Decision Criteria for Determination of Pyrogenicity**

414 The *t*-test is used to compare the data of a test sample against the data of the EC (0.5 EU/mL)
415 that is performed in parallel. If this test results in a significant *p*-value (i.e., smaller than 1%),
416 then the sample is considered to be non-pyrogenic, and as pyrogenic if otherwise (Hoffmann
417 et al., 2005), as long as the assay acceptability criteria in **Section 8.0** has been met.

418 **10.0 STUDY REPORT**

419 The test report should include the following information:

420 *Test Substances and Control Substances*

- 421 • Name of test substance
- 422 • Purity and composition of the substance or preparation
- 423 • Physicochemical properties (e.g., physical state, water solubility)

- 424 • Treatment of the test/control substances prior to testing (e.g., vortexing,
425 sonication, warming, resuspension solvent)

426 *Justification of the In Vitro Test Method and Protocol Used*

427 *Test Method Integrity*

- 428 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
429 test method over time
- 430 • If the test method employs proprietary components, documentation on the
431 procedure used to ensure their integrity from “lot-to-lot” and over time
- 432 • The procedures that the user may employ to verify the integrity of the
433 proprietary components

434 *Criteria for an Acceptable Test*

- 435 • Acceptable concurrent positive control ranges based on historical data
- 436 • Acceptable negative control data

437 *Test Conditions*

- 438 • Cell system used
- 439 • Calibration information for the spectrophotometer used to read the ELISA
- 440 • Details of test procedure
- 441 • Description of any modifications of the test procedure
- 442 • Reference to historical data of the model
- 443 • Description of evaluation criteria used

444 *Results*

- 445 • Tabulation of data from individual test samples

446 *Description of Other Effects Observed*

447 *Discussion of the Results*

448 *Conclusion*

449 *A Quality Assurance Statement for Good Laboratory Practice (GLP)-Compliant Studies*

- 450 • This statement should indicate all inspections made during the study and the
451 dates any results were reported to the Study Director. This statement should
452 also confirm that the final report reflects the raw data.

453 If GLP-compliant studies are performed, then additional reporting requirements provided in
454 the relevant guidelines (e.g., OECD 1998; EPA 2003a, 2003b; FDA 2003) should be
455 followed.

456

456

457 **11.0 REFERENCES**

458 EPA. 2003a. Good Laboratory Practice Standards. Toxic Substances Control Act. 40 CFR
459 792.

460 EPA. 2003b. Good Laboratory Practice Standards. Federal Insecticide, Fungicide, and
461 Rodenticide Act. 40 CFR 160.

462 FDA. 2003. Good Laboratory Practices for Nonclinical Laboratory Studies. 21 CFR 58.

463 Fennrich S, Wendel A, Hartung T. 1999a. New applications of the human whole blood
464 pyrogen assay (PyroCheck). ALTEX 16:146-149.

465 Fennrich S, Fischer M, Hartung T, Lexa P, Montag-Lessing T, Sonntag H-G, Weigand M,
466 Wendel A. 1999b. Detection of endotoxins and other pyrogens using human whole blood.
467 Dev Biol Stand. 101:131-139.

468 Hartung T, Wendel A. 1999. Detection of pyrogens using human whole blood. *In Vitro*
469 Toxicol. 9(4):353-359.

470 Hoffmann S, Peterbauer A, Schindler S, Fennrich S, Poole S, Mistry Y, Montag-Lessing T,
471 Spreitzer I, Löschner B, van Aalderen M, Bos R, Gommer M, Nibbeling R, Werner-
472 Felmayer G, Loitzl P, Jungi T, Brcic M, Brügger P, Frey E, Bowe G, Casado J, Coecke S, de
473 Lange J, Mogster B, Næss LM, Aaberge IS, Wendel A, Hartung T. 2005. International
474 validation of novel pyrogen tests based on human monocytoid cells. J Immunol Methods
475 298:161-173.

476 Organization for Economic Cooperation and Development (OECD). 1998. OECD Series on
477 Principle of Good Laboratory Practice and Compliance Monitoring. No. 1. OECD Principles
478 of Good Laboratory Practice (as revised in 1997). Organisation for Economic Co-operation
479 and Development (OECD), ENV/MC/CHEM (98)17. Paris: OECD.

480 Schindler S, Asmus S, von Aulock S, Wendel A, Hartung T, Fennrich S. 2004.

481 Cryopreservation of human whole blood for pyrogenicity testing. J Immunol Methods
482 294:89-100.

483 Schindler S, Spreitzer I, Löschner B, Hoffmann S, Hennes K, Halder M, Brügger P, Frey E,
484 Hartung T, Montag T. 2006. International validation of pyrogen tests based on cryopreserved
485 human primary blood cells. *J Immunol Methods* 316:42-51.

486

486

487 **12.0 TERMINOLOGY AND FORMULA**488 **12.1 Assay Sensitivity (λ)**

489 For an *in vitro* cell-based assay, the variable λ is defined as the lowest statistically significant
490 point on the standard endotoxin concentration-response curve and represents the relative
491 sensitivity of the test method for the detection of endotoxin (i.e., level of detection).

492 **12.2 Endotoxin Control (EC)**

493 The EC is incubated with Cryo WB and serves as the positive control for the experiment. The
494 results should be compared to historical values to insure that it provides a known level of
495 cytokine release relative to the NSC.

496 **12.3 Endotoxin Limit Concentration (ELC)**

497 The ELC is the maximum allowable concentration of endotoxin for a particular product and
498 is expressed in EU per volume (mL) or weight (mg). The ELC is defined by the FDA or
499 specified in the USP². It is calculated as the product of K/M, where:

500 K is the threshold pyrogen dose for parenteral use in rabbits or humans (5.0 EU/kg). At an
501 injection volume of 10 mL/kg, K is equal to 0.5 EU/mL.

502 M is the larger of the rabbit dose or the maximum human dose administered in one hour as
503 defined below and varies with test substance³.

504 **12.4 Maximum Valid Dilution (MVD)**

505 The MVD is the maximum dilution of a test substance that can be tolerated in a test system
506 without exceeding the ELC, if the test substance must be diluted as a result of assay
507 interference. Dilutions beyond the MVD would not be valid for endotoxin detection in the

² ELC values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, Pharmacopoeial guidelines, FDA publications).

³ Values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, pharmacopoeial guidelines, FDA publications).

508 test system. Calculation of the MVD is dependent on whether or not the ELC for a test
509 substance is published. When the ELC is known, the MVD is:

$$510 \quad \text{MVD} = (\text{ELC} \times \text{Product Potency [PP]})/\lambda$$

511 As an example, for "Cyclophosphamide Injection," the ELC is 0.17 EU/mg, PP is 20
512 mg/mL, and the assay sensitivity is 0.1 EU/mL. The calculated MVD would be 34. The test
513 substance can be diluted no more than 1:34 prior to testing.

514 If the ELC is not known, the MVD is:

$$515 \quad \text{MVD} = \text{PP}/\text{Minimum Valid Concentration (MVC)}$$

$$516 \quad \text{where, MVC} = (\lambda \times \text{M})/\text{K}$$

517 where, M is the maximum human dose

518 As an example, for "Cylophosphamide Injection," the PP is 20 mg/mL, M is 30
519 mg/kg, and assay sensitivity is 0.1 EU/mL. The calculated MVC is 0.6 mg/mL and the MVD
520 is 33.3. The test substance can be diluted no more than 1:33 in the assay prior to testing.

521 **12.5 Negative Product Control (NPC)**

522 The NPC is a test sample to which PFS is added. The NPC is the baseline for determination
523 of cytokine release relative to the endotoxin-spiked PPC.

524 **12.6 Negative Saline Control (NSC)**

525 The NSC is Cryo WB (in RPMI) incubated with PFS (used for dilution of test substance) and
526 is used as the blank.

527 **12.7 Parenteral Threshold Pyrogen Dose (K)**

528 The value K represents the threshold pyrogen dose for parenteral products for rabbits and
529 humans. Based on experimental data, K is fixed at 5.0 EU/kg. For intrathecal products, K is
530 0.2 EU/kg.

531 **12.8 Positive Product Control (PPC)**

532 The PPC is a test substance spiked with the control standard endotoxin (i.e., 0.5 EU/mL or an
533 amount of endotoxin equal to that which produces $\frac{1}{2}$ the maximal increase in OD from the

534 endotoxin standard curve) to insure that the test system is capable of endotoxin detection in
535 the product as diluted in the assay.

536 **12.9 Product Potency (PP)**

537 The concentration for a test substance is the PP typically expressed as $\mu\text{g/mL}$ or mg/mL .

538 **12.10 Rabbit Pyrogen Test Dose or Maximum Human Dose (M)**

539 The variable M represents the rabbit test dose or the maximum human dose in 1 hr. The
540 variable M is expressed in mg/kg and varies with the test substance. For
541 radiopharmaceuticals, M should be adjusted to account for product activity (radioactive
542 decay) at time administration. An average human standard weight of 70 kg is used for the
543 calculation. If a pediatric dose should be used and it is higher than the adult dose, then it
544 should be used in the formula.

545

546