

## **Appendix B3**

### **Draft ICCVAM Recommended Protocol for Future Studies Using the Peripheral Blood Mononuclear Cell (PBMC)/Interleukin-6 (IL-6) Test Method**

***[This Page Intentionally Left Blank]***

1 **Draft ICCVAM Recommended Protocol for Future Studies Using the Peripheral Blood**  
2 **Mononuclear Cell (PBMC)/Interleukin-6 (IL-6) Test Method**

3

4

**PREFACE**

5

6 This proposed protocol for the detection of pyrogenicity is based on information obtained  
7 from 1) The European Centre for the Validation of Alternative Methods (ECVAM)  
8 PBMC/IL-6 Background Review Document (BRD) presented in Appendix A of the draft  
9 Interagency Coordinating Committee on the Validation of Alternative Methods (ICCVAM)  
10 Pyrogenicity Test Method BRD, which includes ECVAM Standard Operating Procedures  
11 (SOPs) for the PBMC/IL-6 test method, and 2) Information provided to the National  
12 Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative  
13 Toxicological Methods (NICEATM) by Dr. Thomas Hartung, Head of ECVAM. The  
14 ECVAM SOPs are based on various test methods that use human PBMCs to detect cytokine  
15 production as a measure of pyrogenicity (Bleeker et al., 1994; Dinarello et al., 1984; Poole et  
16 al., 2003). A table of comparison between the draft ICCVAM recommended protocol and the  
17 ECVAM SOPs is provided in **Table 1**. Future studies using the PBMC/IL-6 test method may  
18 include further characterization of the usefulness or limitations of the assay for regulatory  
19 decision-making. Users should be aware that the proposed test method protocol might be  
20 revised based on additional optimization and/or validation studies. ICCVAM recommends  
21 that test method users routinely consult the ICCVAM/NICEATM website  
22 (<http://iccvam.niehs.nih.gov/>) to ensure that the most current test method protocol is used.

23

23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38

***[This Page Intentionally Left Blank]***

39 **Table 1 Comparison of Draft Recommended Test Method Protocol with the ECVAM SOP for the Peripheral Blood**  
 40 **Mononuclear Cell (PBMC)/Interleukin-6 (IL-6) Pyrogen Test**

Protocol Component	Draft ICCVAM Protocol	ECVAM SOP/ ECVAM Catch-Up Validation SOP
<b>Test Substance</b>	Test neat or at minimal dilution that produces no interference	Test at MVD during validation - otherwise similar to ICCVAM protocol
<b>Incubation Plate (number of control or test groups at n=4 each)</b>	NSC (1)	Same as ICCVAM protocol
	EC (5)	Same as ICCVAM protocol
	TS (14)	Same as ICCVAM protocol
<b>ELISA Plate</b>	Includes seven point IL-6 SC and blank in duplicate	Same as ICCVAM protocol
<b>Decision Criteria for Interference</b>	$0.5 \times \text{Median OD}_{450}^1$ of 1 EU/mL EC <2x Median $\text{OD}_{450}$ of 1 EU/mL EC	$0.5 \times \text{Median OD}_{450}$ of 0.25 EU/mL <sup>3</sup> EC <2x Median $\text{OD}_{450}$ of 0.25 EU/mL EC
<b>Assay Acceptability Criteria</b>	Mean $\text{OD}_{450}$ of PPC is 50% to 200% of 0.5 EU/mL EC	Mean $\text{OD}_{450}$ of PPC is 50% to 200% of 0.25 EU/mL <sup>2</sup> EC
	Mean $\text{OD}_{450}$ of NSC $\leq 0.15$	Same as ICCVAM protocol
	Quadratic function of IL-6 SC $r^2 \geq 0.95$	Same as ICCVAM protocol
	EC SC produces $\text{OD}_{450}$ values that ascend in a sigmoidal concentration response	Same as ICCVAM protocol
	High responder (i.e., > 200 pg/mL IL-6) or low responder (i.e., Mean $\text{OD}_{450}$ for 1 EU/mL EC is significantly less than that for 1000 pg/mL IL-6) blood donors may be excluded (limited to 1 of 4 donors). One of four donors may be excluded or the test repeated with four additional donors.	Same as ICCVAM protocol
<b>Decision Criteria for Pyrogenicity</b>	Mean corrected $\text{OD}_{450}$ of TS > Mean corrected $\text{OD}_{450}$ of 0.5 EU/mL EC	EC SC data transformed to 4-parameter logistical model by an in-house program and the linear mean square is calculated. TS pyrogen content is compared with the ELC <sup>3</sup> using confidence limits for significance. The preparation being examined must pass the test with blood from three separate donors. Dixon's test is used to reject outliers.

Protocol Component	Draft ICCVAM Protocol	ECVAM SOP/ ECVAM Catch-Up Validation SOP
	<p><i>Decision Level 1</i>                      Passes if all are negative                      Fails if 2 or more are positive                      If one is positive, repeat with 4 more donors</p> <p><i>Decision Level 2</i>                      Passes if 1/8 are positive                      Fails if 2 or more of 8 are positive</p>	<p>Same as ICCVAM protocol</p>
	<p>Not included</p>	<p>Limit test is run to determine whether or not a TS after correction and dilution contains &lt; 0.5 EU/mL of endotoxin</p>

41 Abbreviations: EC = Endotoxin control; ELC = Endotoxin Limit Concentration; MVD = Maximum Valid Dilution; NSC = Normal saline control; PPC =  
 42 Positive Product Control; SC = Standard curve; TS = Test substance

43 <sup>1</sup> Median or mean OD<sub>450</sub> values are corrected (i.e., reference filter reading, if applicable, and NSC are subtracted).

44 <sup>2</sup> Value of 0.25 EU/mL is arbitrary - should correspond to a concentration in the middle of the SC for the EC.

45 <sup>3</sup> Where unknown, the ELC is calculated (e.g., Based on a rabbit sensitivity of 5 EU/kg, for a product injected at 10 mL/kg, the detection limit is 5 EU/10 mL/kg  
 46 or 0.5 EU/mL/kg, or an ELC of 0.5 EU/mL).

47

48

49

50

51

52

53

54

55

**56 1.0 PURPOSE AND APPLICABILITY**

57 The purpose of this protocol is to describe the procedures used to evaluate the presence of a  
58 pyrogen (i.e., Gram-negative endotoxin) in a test substance. The presence of Gram-negative  
59 endotoxin is detected by its ability to induce cytokine IL-6 release from human PBMCs. The  
60 quantity of IL-6 released is obtained using an enzyme-linked immunosorbent assay (ELISA)  
61 that includes monoclonal or polyclonal antibodies specific for IL-6. Release of this cytokine  
62 is measured by incubation of PBMCs with test substances or controls (i.e., positive,  
63 negative). The amount of pyrogen present is determined by comparing the values of  
64 endotoxin equivalents produced by PBMCs exposed to the test substance to those exposed to  
65 an internationally-harmonized Reference Standard Endotoxin (RSE)<sup>1</sup> or an equivalent  
66 standard expressed in Endotoxin Units (EU)/mL.

67 The focus of this protocol is on the use of the PBMC/IL-6 test method specifically for the  
68 detection of Gram-negative endotoxin in parenteral pharmaceuticals. The relevance and  
69 reliability for non-endotoxin pyrogens (e.g., lipoteichoic acid) has not been demonstrated in a  
70 formal validation study.

**71 2.0 SAFETY AND OPERATING PRECAUTIONS**

72 All procedures for procurement of eligible blood donors and blood donations should follow  
73 the regulations and procedures set forth by institutional guidelines for utilization of human  
74 substances, which include but are not limited to blood, tissues, and tissue fluids. Standard  
75 laboratory precautions are recommended including the use of laboratory coats, eye  
76 protection, and gloves. If necessary, additional precautions required for specific study  
77 substances or hazardous chemicals will be identified in the Material Safety Data Sheet  
78 (MSDS).

---

<sup>1</sup> RSEs are internationally-harmonized reference standards (e.g., WHO *E. coli* Lipopolysaccharide [LPS] 94/580 [0113:H10:K-]; USP RSE Lot G3E069; FDA Lot EC-6). Equivalent endotoxins include commercially available *E. coli*-derived LPS Control Standard Endotoxin (CSE) or other *E. coli* LPS preparations that have been calibrated with an appropriate RSE.

79 The stop solution used in the ELISA kit is acidic and corrosive and should be handled with  
80 the proper personal protective devices. If this reagent comes into contact with skin or eyes,  
81 wash thoroughly with water. Seek medical attention, if necessary.

82 Tetramethylbenzidine (TMB) solution contains a hydrogen peroxide substrate and 3, 3', 5,  
83 5'- TMB. This reagent is a strong oxidizing agent and a suspected mutagen. Appropriate  
84 personal protection should be used to prevent bodily contact.

85 Bacterial endotoxin is a toxic agent (i.e., can induce sepsis, shock, vascular damage,  
86 antigenic response) and should be handled with care. Skin cuts should be covered and  
87 appropriate personal protective devices should be worn. In case of contact with endotoxin,  
88 immediately flush eyes or skin with water for at least 15 min. If inhaled, remove the affected  
89 individual from the area and provide oxygen and/or artificial respiration as needed. Skin  
90 absorption, ingestion, or inhalation may produce fever, headache, and hypotension.

### 91 **3.0 MATERIALS, EQUIPMENT, AND SUPPLIES**

#### 92 **3.1 Source of Cells**

93 PBMCs from fresh WB are the source of cells for cytokine production in the PBMC/IL-6 test  
94 method (Poole et al., 2003). WB is obtained from healthy human volunteers who have  
95 provided their consent according to established institutional guidelines. Volunteers are  
96 expected not to have taken any drugs (e.g., prescription drugs, recreational drugs, herbal  
97 drugs) and to have been free from illness for at least two weeks prior to donation.

98 PBMCs are isolated from WB using density gradient centrifugation (refer to **Section 6.1.1.1**).  
99 The isolated PBMC suspension may be used directly in the PBMC/IL-6 test assay (**Section**  
100 **6.1.3**) or frozen for later use (**Section 6.1.1.3**).

#### 101 **3.2 Equipment and Supplies**

102 For all steps in the protocol, excluding the ELISA procedure, the materials that will be in  
103 close contact with samples and/or blood cells (e.g., pipet tips, containers, solutions) should be  
104 sterile and free from detectable pyrogens.



105 3.2.1 Preparation of PBMCs106 3.2.1.1 *Equipment*

- 107 • Centrifuge
- 108 • Hood; Bio-safety, laminar flow (recommended)
- 109 • Incubator; cell culture (37±1°C + 5% CO<sub>2</sub>)
- 110 • Lymphoprep™
- 111 • Pipetter; multichannel (8- or 12-channel)
- 112 • Pipettors; single-channel adjustable (20, 200, and 1000 µL)
- 113 • Repeating pipetter
- 114 • Vortex mixer

115 3.2.1.2 *Consumables*

- 116 • Centrifuge tubes; nonpyrogenic, polystyrene (15 and 50 mL)
- 117 • Combitips; repeating pipetter (2.5 and 5.0 mL)
- 118 • Cryotubes; screw-cap, 2 mL
- 119 • Filters; sterile, 0.22 µm
- 120 • Needle set; Sarstedt multifly, pyrogen-free, 19 mm, 21 gauge for S-Monovette
- 121 • Phosphate buffered saline (PBS); sterile
- 122 • Pipets; serological, sterile (5, 10, and 25 mL)
- 123 • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 124 • Reaction tubes; polystyrene (1.5 mL)
- 125 • Reservoirs; fluid
- 126 • RPMI-1640 cell culture medium; supplemented with the following reagents to
- 127 yield RPMI-Complete (RPMI-C)
  - 128 ○ FCS; heat-inactivated (5 mL or a 1% final concentration)
  - 129 ○ L-Glutamine; 2 mM

- 130                   ○ Penicillin/streptomycin (10,000 IU/mL penicillin, 10 mg/mL
- 131                   streptomycin)
- 132                   • Syringes; sterile (100 µL and 30 mL)
- 133                   • Tips; pipetter, sterile, pyrogen-free (20, 200, and 1000 µL)
- 134                   • Tubes; Sarstedt S-Monovette, 7.5 mL, heparinized for blood collection

### 135 3.2.2 ELISA

#### 136 3.2.2.1 *Equipment*

- 137                   • Microplate mixer
- 138                   • Microplate reader (450 nm with an optional reference filter in the range of
- 139                   540-590 nm)
- 140                   • Microplate washer (optional)
- 141                   • Multichannel pipetter

#### 142 3.2.2.2 *Consumables*

- 143                   • Container; storage, plastic
- 144                   • Deionized water; nonsterile
- 145                   • Plates; microtiter, nonpyrogenic, 96-well, polystyrene, tissue culture
- 146                   • Pyrogen-free water (PFW)
- 147                   • Reservoirs; fluid
- 148                   • Tips; pipetter, nonsterile
- 149                   • Tubes; polystyrene (12 mL)

#### 150 3.2.2.3 *ELISA Kit*

151 An ELISA that measures IL-6 release from PBMCs is used. A variety of IL-6 ELISA kits are  
152 commercially available and the IL-6 ELISA procedure outlined in this protocol is intended to  
153 serve as an example for using an ELISA kit. If the user prefers to prepare an in-house  
154 ELISA, then additional reagents would be required. The IL-6 ELISA should be calibrated  
155 using an IL-6 international reference standard (e.g., WHO 89/548) prior to use. The IL-6

156 cytokine assay kits do not provide the RSE or endotoxin equivalent; therefore, this reagent  
157 must be purchased separately. Results obtained using these products are subject to the assay  
158 acceptability and decision criteria described in **Sections 8.0** and **9.0**. IL-6 ELISA kit  
159 components may include the following:

- 160 • ELISA plates coated with anti-human IL-6 capture antibody; monoclonal or  
161 polyclonal
- 162 • Buffered wash solution
- 163 • Dilution buffer
- 164 • Enzyme-labeled detection antibody
- 165 • Human IL-6 reference standard
- 166 • Pyrogen-free saline (PFS)
- 167 • Stop solution
- 168 • TMB/substrate solution

### 169 **3.3 Chemicals**

- 170 • Endotoxin (e.g., WHO *E. coli* LPS 2nd International Standard 94/580; USP  
171 RSE *E. coli* LPS Lot G3E069; USP RSE *E. Coli* Lot G; FDA *E. coli* Lot EC6)

### 172 **3.4 Solutions**

- 173 • RPMI-C cell culture medium

## 174 **4.0 ASSAY PREPARATION**

175 All test substances, endotoxin, and endotoxin-spiked solutions should be stored at 4°C.

### 176 **4.1 Endotoxin Standard Curve**

177 An internationally harmonized RSE or equivalent is used to generate the endotoxin standard  
178 curve. The use of any other *E. coli* LPS requires calibration against a RSE using the  
179 PBMC/IL-6 test method.

180 A standard endotoxin curve consisting of a Normal Saline Control (NSC) and five RSE  
 181 concentrations (0.063, 0.125, 0.25, 0.50, and 1.0 EU/mL) are included in the incubation step  
 182 (refer to **Table 4-1**) and then transferred to the ELISA plate. To prepare the endotoxin  
 183 standard curve, first obtain a 2000 EU/mL stock solution by addition of PFW to the  
 184 lyophilized content of the stock vial by following the instructions provided by the  
 185 manufacturer (e.g., for a vial containing 10,000 EU, 5 mL of PFW is added). To reconstitute  
 186 the endotoxin, the stock vial should be vortexed vigorously for at least 30 min or sonicated in  
 187 a bath sonicator for 5 min. The stock solution is stable for 14 days when stored at 2 to 8°C.  
 188 An endotoxin standard curve is prepared by making serial dilutions of the stock solution in  
 189 PFS as described in **Table 4-1**.

190 **Table 4-1 Preparation of Endotoxin Standard Curve**

Stock Endotoxin EU/mL	µL of Stock Endotoxin	µL of PFS	Endotoxin Concentration EU/mL
2000 <sup>1,2</sup>	20	1980	20 <sup>3</sup>
20	100	1900	1.0
1.0	500	500	0.50
0.50	500	500	0.25
0.25	500	500	0.125
0.125	500	500	0.063
0	0	1000	0

191 Abbreviations: EU = Endotoxin units; PFS = Pyrogen-free saline

192 Each stock tube should be vortexed vigorously prior to its use to make the subsequent dilution.

193 <sup>1</sup> A 2000 EU/mL stock solution of endotoxin is prepared according to the manufacturer's instructions.

194 <sup>2</sup> The stock solution of USP RSE may be stored in aliquots and kept at -20°C for up to 6 months. Do not store  
 195 the endotoxin at -80°C.

196 <sup>3</sup> This concentration is not used in the assay.

197

## 198 **4.2 Test Substances**

199 Liquid test substances should be tested neat or, if interference is detected (see **Section 4.2.1**),  
 200 diluted in PFS. Solid test substances should be prepared as solutions in PFS or, if insoluble in  
 201 saline, dissolved in dimethylsulfoxide (DMSO) then diluted up to 0.5% (v/v) with PFS  
 202 provided that this concentration does not interfere with the assay. The test substances should  
 203 be vortexed vigorously for at least 30 min or sonicated in a bath sonicator for 5 min.

#### 204 4.2.1 Interference Testing

205 Interference testing must be carried out on any test sample for which no interference  
206 information is available. The purpose of the interference test is to determine the lowest  
207 dilution (i.e., highest concentration) of a test substance from which an endotoxin spike can be  
208 detected (i.e., based on the decision criteria described in **Section 4.2.1.2**). However, to ensure  
209 a valid test, a test substance should not be diluted beyond its Maximum Valid Dilution  
210 (MVD).

211 For many marketed products, values for the MVD and the Endotoxin Limit Concentration  
212 (ELC) are published in the U.S. Pharmacopeia, the European Pharmacopoeia, and/or Food  
213 and Drug Administration (FDA) guidelines. However if one or both of these values are not  
214 available, then calculation of the MVD is dependent on the ELC (see **Section 12.3**). If  
215 unknown, the ELC can be approximated by dividing the maximum hourly dose of the  
216 product by the hourly dose received per patient. For example, if a product is used at an  
217 hourly dose of 100 mg per patient, then the ELC would be 350 EU/100 mg or 3.5 EU/mg.

##### 218 4.2.1.1 *Reference Endotoxin for Spiking Test Substances*

219 The WHO-LPS 94/580 [*E. coli* O113:H10:K-] or equivalent international RSE (e.g., USP  
220 G3E069, FDA EC-6) is recommended for preparation of the endotoxin control (EC). If a  
221 different *E. coli* LPS is used and the potency relative to the RSE is not provided, then each  
222 lot must be calibrated against the RSE in the PBMC/IL-6 test method. For interference  
223 testing, an endotoxin standard curve (see **Section 4.1**) should be included on each plate.

##### 224 4.2.1.2 *Spiking Test Substances with Endotoxin*

225 For interference testing, non-spiked and endotoxin-spiked test substances are prepared in  
226 microplate wells (n=4 replicates) and an *in vitro* pyrogen test is performed. Either RPMI-C or  
227 a fixed concentration (a concentration selected from the middle of the EC standard curve) of  
228 the RSE (i.e., 0.25 EU/mL) in RPMI-C is added to the test substance in serial two-fold  
229 dilutions. An illustrative example of endotoxin spiking solutions is shown in **Table 4-2**. For  
230 non-spiked solutions, 150  $\mu$ l of RPMI-C is added to a well followed by 50  $\mu$ l of the test  
231 substance (neat or at serial dilution) and 50  $\mu$ L of PBMCs and the well contents are mixed.  
232 Endotoxin-spiked solutions are prepared by adding 100  $\mu$ L of RPMI-C to each well followed

233 by 50  $\mu\text{L}$  of the test substance (neat or at serial dilution) and 50  $\mu\text{L}$  of an endotoxin spike  
 234 solution (0.25 EU/mL). Finally, 50  $\mu\text{L}$  of PBMCs are added and the well contents are mixed  
 235 (see example presented in **Table 4-2**).

236 **Table 4-2 Preparation of Endotoxin-Spiked and Non-spiked Solutions for**  
 237 **Determination of Test Substance Interference in the Incubation and**  
 238 **ELISA Test Systems**

Sample Addition	Spiked	Non-spiked
	$\mu\text{L}/\text{well}^1$	
RPMI-C	100	150
Endotoxin spike solution <sup>2</sup>	50	0
Test substance (neat and each serial dilution)	50	50
PBMCs <sup>3</sup>	50	50
Total <sup>4</sup>	250	250

239 Abbreviations: PBMC = Peripheral blood mononuclear cells

240 <sup>1</sup> n=4 replicates each

241 <sup>2</sup> Endotoxin concentration is 0.25 EU/mL in RPMI-C.

242 <sup>3</sup> PBMCs are resuspended in RPMI-C ( $1 \times 10^6$  cells/mL).

243 <sup>4</sup> A total volume of 250  $\mu\text{L}$  per well is used for the incubation.

244

245 The lowest dilution of the test substance that yields an endotoxin spike recovery of 50% to  
 246 200% in the pyrogen test is determined. The optical density (OD) values of the endotoxin-  
 247 spiked and non-spiked test substances are calibrated against the endotoxin calibration curve.  
 248 The resulting EU value of the non-spiked test substance is subtracted from the corresponding  
 249 EU value of the endotoxin-spiked test substance at each dilution. The % recovery for each  
 250 sample dilution is then determined from the endotoxin spike solution concentration set to  
 251 100%. For example, consider the following interference test results in **Table 4-3**:

252 **Table 4-3 Example of Interference Data Used to Determine Sample Dilution for**  
 253 **Assay**

Sample Dilution	% Recovery of Endotoxin Control
None	25
1:2	49
1:4	90
1:8	110

254

255 Based on these results, the dilution of the test substance used in the *in vitro* pyrogen test  
 256 would be 1:4 (i.e., the lowest dilution between 50% and 200% of the 0.25 EU/mL EC).

#### 257 4.2.2 Interference with ELISA

258 If the data obtained from the experiment in **Section 4.2.1** suggests the presence of  
259 interference, then a subsequent experiment similar to that described in **Section 4.2.1** would  
260 need to be performed to confirm that the test substance(s) does not directly interfere with the  
261 ELISA. For this experiment, an ELISA would be performed in the absence of PBMCs.

### 262 **5.0 CONTROLS**

#### 263 **5.1 Negative Control**

264 A negative control (e.g., RPMI-C) is included in each experiment in order to detect  
265 nonspecific changes in the test system, as well as to provide a baseline for the assay  
266 endpoints.

#### 267 **5.2 Solvent Control**

268 Solvent controls are recommended to demonstrate that the solvent is not interfering with the  
269 test system when solvents other than PFS are used to dissolve test substances.

#### 270 **5.3 Positive Control**

271 An internationally standardized EC (e.g., WHO 94/580; USP G3E069; 0.5 EU/mL) is  
272 included in each experiment to verify that an appropriate response is induced.

#### 273 **5.4 Benchmark Control**

274 Benchmark controls may be used to demonstrate that the test method is functioning properly,  
275 or to evaluate the relative pyrogenic potential of chemicals (e.g., parenteral pharmaceuticals,  
276 medical device eluates) of a specific class or a specific range of responses, or for evaluating  
277 the relative pyrogenic potential of a test substance. Appropriate benchmark controls should  
278 have the following properties:

- 279 • consistent and reliable source(s) for the chemicals (e.g., parenteral  
280 pharmaceuticals, medical device eluates)
- 281 • structural and functional similarities to the class of the substance being tested
- 282 • known physical/chemical characteristics

- 283                   • supporting data on known effects in animal models
- 284                   • known potency in the range of response

## 285   **5.5           Positive Product Control (PPC)**

286   The PPC is a test substance diluted to a level that does not interfere with the test method and  
287   does not exceed the MVD. The PPC is obtained by spiking a test substance with a known EC  
288   (e.g., 0.25 EU/mL) and demonstrating that 50% to 200% of the EC is recovered.

## 289   **5.6           Negative Product Control (NPC)**

290   The NPC is the test substance diluted to the MVD and then spiked with PFS. It is the  
291   negative control for the PPC.

## 292   **6.0           EXPERIMENTAL DESIGN**

### 293   **6.1           Incubation with Test Samples and Measurement of IL-6 Release**

#### 294   6.1.1       Collection of Human Blood

295   WB is obtained from healthy human volunteers who have provided their consent according  
296   to established institutional guidelines. Volunteers are expected not to have taken any drugs  
297   and to have been free from illness for at least two weeks prior to donation. The criteria for  
298   rejection of data from donors that are low responders or that are suspect due to veracity of  
299   health information is addressed in **Section 8.0**.

300   Fresh WB is drawn by venipuncture using a multily needle set and collected in heparinized  
301   tubes (15 IU/mL lithium heparin). All components of the blood collection system (i.e.,  
302   syringes, tubes, connecting lines) must be sterile and pyrogen-free. Blood should be stored at  
303   room temperature (RT) and must be used within 4 hr. Prior to use in the assay, the collection  
304   tubes should be gently inverted 1 to 2 times. **Do Not Vortex**.

#### 305   6.1.1.1    *Isolation of PBMCs from WB*

306   PBMCs must be isolated within two hours of WB collection using Lymphoprep™. The  
307   isolation procedure described below is a modification of the manufacturer's instructions as  
308   outlined in the ECVAM SOPs for the PBMC/IL-6 test method.



309 To each 15 mL of heparinized WB, add 15 mL of PBS and a sterile 10 mL pipet into the  
310 tube. Add 20 mL of Lymphoprep™ to each of the two tubes, through the inserted pipet, to  
311 form a lower, denser, layer. The tubes are then centrifuged at 340 x g for 45 min at room RT.  
312 After centrifugation, a white band of PBMCs should be visible at approximately the 25 mL  
313 graduation mark on the tube. Carefully remove the uppermost 9 mL from each tube and  
314 discard. If cryopreservation of PBMCs is to be performed (see **Section 6.1.1.3**), then transfer  
315 the layer to a new 50 mL tube for preparing a cryoprotective solution. The remaining  
316 supernatant above the PBMC band should be aspirated and discarded. Using a 10 mL pipet,  
317 transfer the PBMC layer to a new 50 mL centrifuge tube.

#### 318 6.1.1.2 *Washing PBMCs*

319 The PBMCs are resuspended in a total volume of 50 mL of PBS and centrifuged at 340 x g  
320 for 15 min. The supernatant is poured off and the cellular sediment resuspended in 10 mL of  
321 PBS by pipetting up and down several times with a serological pipet. Adjust the total volume  
322 in each tube to 50 mL with PBS and centrifuge at 340 x g for 10 min for a second time. After  
323 centrifugation, the PBMCs are resuspended with 15 mL of RPMI-C and the two aliquots of  
324 PBMCs are pooled into a single tube. Cell counts (expressed per mL) are determined using a  
325 hemacytometer. If the cell count is above  $1.2 \times 10^6/\text{mL}$ , the cell suspension should be diluted  
326 to  $1 \times 10^6$  cell/mL in RPMI-C. This suspension must be used in the PBMC/IL-6 assay within  
327 4 hr from the time of WB collection.

#### 328 6.1.1.3 *Procedure for Cryopreservation and Thawing of PBMCs*

329 To freeze the PBMCs, prepare a cryoprotective solution by adding 2 mL of pyrogen-free  
330 DMSO to the 18 mL of supernatant collected in the centrifugation procedure outlined in  
331 **Section 6.1.1.1**. Cool the cryoprotective solution to between 2 and 8°C. Centrifuge the  
332 isolated PBMCs as instructed in **Section 6.1.1.2** and then add 6 mL of the chilled  
333 cryoprotective solution to the cell sediment. Pool the cell suspensions from the same donor  
334 and transfer 1.0 mL aliquots to appropriately labeled screw-cap cryotubes. The cryotubes are  
335 placed in a styrofoam box for thermal insulation and slowly frozen to -80°C. After 72 hr, the  
336 tubes can be transferred to liquid nitrogen for prolonged storage.

337 To thaw the cryopreserved PBMCs, take two tubes (each at 1.0 mL) from a single donor and  
338 submerge in a water bath at  $37 \pm 1^\circ\text{C}$ . After thawing, the cell suspensions are pooled in a

339 single 50 mL centrifuge tube and RPMI-C is added to give a total volume of 40 mL. The  
 340 PMBCs are centrifuged at 340 x g for 10 min, the supernatant removed, and the cells  
 341 resuspended in 10 mL of RPMI-C. Viability may be tested and should yield greater than 95%  
 342 cell survival.

### 343 6.1.2 Incubation Plate

344 Test substances are prepared at a level of dilution that did not show interference with the test  
 345 system or for which it is known that interference does not occur. Each incubation plate can  
 346 accommodate an endotoxin standard curve, a NSC, and 14 test samples (see **Table 6-1**).

347 **Table 6-1 Overview of Incubation Plate Preparation in the PBMC/IL-6 Test**  
 348 **Method**

Number of Wells	Sample	RPMI-C	EC	Test Sample	PBMCs	Mix the samples; incubate overnight at 37±1°C in a humidified atmosphere with 5% CO <sub>2</sub> .	Mix the samples; immediately transfer to an ELISA plate <sup>3</sup> and run ELISA or store plate at -20°C or -80°C.
		µL					
20 <sup>1</sup>	EC	100	50	0	100		
4	NSC	150	0	0	100		
56 <sup>2</sup>	Test samples (1-14)	100	0	50	100		

349 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; PBMC = Peripheral blood  
 350 mononuclear cell

351 <sup>1</sup> Five EC concentrations (0.063, 0.125, 0.25, 0.50, and 1.0 EU/mL) in quadruplicate.

352 <sup>2</sup> 14 test samples (n=4 each) per plate.

353 <sup>3</sup> An IL-6 standard curve is prepared in Columns 11 and 12 on the ELISA plate. Therefore, 80 wells are  
 354 available for test samples and controls on the incubation plate.

355

### 356 6.1.3 Incubation Assay for IL-6 Release

357 PBMC samples are prepared in a microtiter plate using a laminar flow hood. All  
 358 consumables and solutions must be sterile and pyrogen-free. Each plate should be labeled  
 359 appropriately with a permanent marker. An overview of the incubation plate preparation is  
 360 shown in **Table 6-1**. The incubation procedure is outlined below:

361 **Step 1.** Refer to the suggested incubation plate template presented in **Table 6-2**.

362 **Step 2.** Using a pipetter, transfer 100 µL of RPMI-C into each well.

363 **Step 3.** Transfer 50  $\mu$ L of test sample into the appropriate wells as indicated in the  
364 template.

365 **Step 4.** Transfer 50  $\mu$ L of the EC (standard curve) and the NSC controls in  
366 quadruplicate into the appropriate wells according to the template.

367 **Step 5.** Transfer 100  $\mu$ L of a well-mixed PBMC suspension into each well and mix  
368 by gently swirling the plate.

369 **Step 6.** Mix the contents of the wells thoroughly by pipetting up and down several  
370 times using a multichannel pipetter, changing the tips between each row in order to  
371 avoid cross-contamination.

372 **Step 7.** Place the covered plate in a tissue culture incubator for 16 to 24 hr at  
373  $37\pm 1^\circ\text{C}$  in a humidified atmosphere containing 5%  $\text{CO}_2$ .

374 **Step 8.** Prior to transferring the test samples onto the ELISA plate, mix the  
375 contents of the wells by pipetting up and down using a multichannel pipetter,  
376 changing the tips between each row in order to avoid cross-contamination.

377 **Table 6-2 Incubation Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	EC <sup>1</sup> 1.0	EC 1.0	EC 1.0	EC 1.0	TS3	TS3	TS3	TS3	TS11	TS11	Void <sup>3</sup>	Void
<b>B</b>	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS4	TS4	TS4	TS4	TS11	TS11	Void	Void
<b>C</b>	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS5	TS5	TS5	TS5	TS12	TS12	Void	Void
<b>D</b>	EC 0.125	EC 0.125	EC 0.125	EC 0.125	TS6	TS6	TS6	TS6	TS12	TS12	Void	Void
<b>E</b>	EC 0.063	EC 0.063	EC 0.063	EC 0.063	TS7	TS7	TS7	TS7	TS13	TS13	Void	Void
<b>F</b>	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	Void	Void
<b>G</b>	TS1 <sup>2</sup>	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	Void	Void
<b>H</b>	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	Void	Void

378 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

379 <sup>1</sup> EC value (e.g., EC 1.0) represents the endotoxin concentration in EU/mL.

380 <sup>2</sup> TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.

381 <sup>3</sup> Columns 11 and 12 are reserved for the IL-6 standard curve on the ELISA plate (see **Table 6-3**).

382

383

384 **6.2 ELISA to Measure IL-6 Release**385 6.2.1 IL-6 Standard Curve

386 An IL-6 standard supplied with the ELISA kit is used. IL-6 standards are typically supplied  
387 in lyophilized form and should be reconstituted according to the manufacturer's instructions.  
388 The stock solution should be diluted in RPMI-C to the following concentrations: 0, 62.5, 125,  
389 250, 500, 1000, 2000, and 4000 pg/mL in volumes of at least 500  $\mu$ L. Each well on the  
390 ELISA plate will receive 50  $\mu$ L of an IL-6 blank or standard.

391 6.2.2 ELISA

392 The manufacturer's instructions provided with the ELISA kit should be followed and a  
393 typical experimental design is outlined below. If the user prefers to prepare an in-house  
394 ELISA, then appropriate modifications and validation of these changes would be necessary.  
395 The ELISA should be carried out at RT and therefore all components must be at RT prior to  
396 use. Do *not* thaw frozen specimens by heating them in a water bath. A suggested ELISA  
397 plate template is shown in **Table 6-3**, which includes a five-point EC standard curve, a NSC,  
398 an eight-point IL-6 standard curve (0 to 4000 pg/mL), and 14 test substances in  
399 quadruplicate. The EC standard curve, the NSC, and the test sample supernatants are  
400 transferred directly from the incubation plate. The IL-6 standard curve is prepared as  
401 described in **Section 6.2.1**. An overview of the ELISA plate preparation is shown in **Table 6-**  
402 **4**.

403 **Step 1.** After pipetting up and down very carefully three times (avoid detachment  
404 of the adherent PBMCs) to mix the supernatant, transfer 50  $\mu$ L from each well of  
405 the Incubation Plate (A1-10; H1-10) to the ELISA plate.

406 **Step 2.** Add 50  $\mu$ L of each IL-6 standard (0 to 4000 pg/mL) into the respective  
407 wells on the ELISA plate.

408 **Step 3.** Add 200  $\mu$ L of the enzyme-labeled detection antibody (neat as supplied, or  
409 diluted, if necessary) to each of the wells.

410           **Step 4.** Cover the microtiter plate(s) with adhesive film and incubate for 2 to 3 hr  
411           at 20 to 25°C.

412           **Step 5.** Decant and wash each well three times with 300 µL Buffered Wash  
413           Solution and then rinse three times with deionized water. Place the plates upside  
414           down and tap to remove water.

415           **Step 6.** Add 200 µL of TMB/Substrate Solution to each well and incubate at RT in  
416           the dark for 15 min. If necessary, decrease the incubation time.

417           **Step 7.** Add 50 µL of Stop Solution to each well.

418           **Step 8.** Tap the plate gently after the addition of Stop Solution to aid in mixing.

419           **Step 9.** Read the OD<sub>450</sub> within 15 min of adding the Stop Solution. Measurement  
420           with a reference wavelength of 540 to 590 nm is recommended.

421

421

422 **Table 6-3 ELISA Plate - Sample and Control Template**

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	EC <sup>1</sup> 1.0	EC 1.0	EC 1.0	EC 1.0	TS3	TS3	TS3	TS3	TS11	TS11	IL-6 <sup>3</sup> 0	IL-6 0
<b>B</b>	EC 0.50	EC 0.50	EC 0.50	EC 0.50	TS4	TS4	TS4	TS4	TS11	TS11	IL-6 62.5	IL-6 62.5
<b>C</b>	EC 0.25	EC 0.25	EC 0.25	EC 0.25	TS5	TS5	TS5	TS5	TS12	TS12	IL-6 125	IL-6 125
<b>D</b>	EC 0.125	EC 0.125	EC 0.125	EC 0.125	TS6	TS6	TS6	TS6	TS12	TS12	IL-6 250	IL-6 250
<b>E</b>	EC 0.063	EC 0.063	EC 0.063	EC 0.063	TS7	TS7	TS7	TS7	TS13	TS13	IL-6 500	IL-6 500
<b>F</b>	NSC	NSC	NSC	NSC	TS8	TS8	TS8	TS8	TS13	TS13	IL-6 1000	IL-6 1000
<b>G</b>	TS1 <sup>2</sup>	TS1	TS1	TS1	TS9	TS9	TS9	TS9	TS14	TS14	IL-6 2000	IL-6 2000
<b>H</b>	TS2	TS2	TS2	TS2	TS10	TS10	TS10	TS10	TS14	TS14	IL-6 4000	IL-6 4000

423 Abbreviations: EC = Endotoxin control; NSC = Negative saline control; TS = Test substance

424 <sup>1</sup> EC value (e.g., EC 1.0) represents the endotoxin concentration in EU/mL.425 <sup>2</sup> TS number (e.g., TS1) represents an arbitrary sequence for individual test substances.426 <sup>3</sup> IL-6 values in columns 11 and 12 are in pg/mL.

427

428

429 **Table 6-4 Overview of ELISA Procedure**

<b>Material transfer from Incubation Plate (µL)</b>	<b>Enzyme-labeled Antibody (µL)</b>	Cover the Incubation Plate and incubate for 2 to 3 hr at 20-25°C.	<b>TMB/Substrate Solution (µL)</b>	Incubate 15 min at RT in dark.	<b>Stop Solution (µL)</b>	Read optical density at 450 nm with a 540-590 nm wavelength reference filter.
50	200		200		50	

430 Abbreviations: RT = Room temperature

431

432

432

433 **7.0 EVALUATION OF TEST METHODS**434 **7.1 OD Measurements**

435 The OD of each well is obtained by reading the samples in a standard microplate  
436 spectrophotometer (i.e., plate reader) using a visible light wavelength of 450 nm (OD<sub>450</sub>) with  
437 a 540 to 590 nm reference filter (recommended). OD<sub>450</sub> values are used to determine assay  
438 acceptability and in the decision criteria for the detection of endotoxin in a test substance (see  
439 **Sections 8.0 and 9.0**).

440 **8.0 CRITERIA FOR AN ACCEPTABLE TEST**

441 Obtain the PPC and the corresponding NPC by interference testing of a test substance in the  
442 presence and absence of a fixed quantity of endotoxin (i.e., 0.25 EU/mL) in quadruplicate.  
443 An EC (five-point standard curve) and a NSC should be included in each experiment. An IL-  
444 6 standard curve should be included in each ELISA as shown in the template presented in  
445 **Table 6-3**. An assay is considered acceptable only if the following minimum criteria are met:

- 446 • The quadratic function of the IL-6 standard curve produces an  $r \geq 0.95$  and the  
447 OD<sub>450</sub> of the blank control is below 0.15.
- 448 • The endotoxin standard curve produces OD<sub>450</sub> values that ascend in a  
449 sigmoidal concentration response.

450 Blood donors are considered to be low responders if their mean OD<sub>450</sub> value obtained for 1  
451 EU/mL EC is significantly below the mean OD<sub>450</sub> value obtained for 1000 pg/mL IL-6.

452 Blood donors who produce a mean OD<sub>450</sub> value for the NSC that is significantly above the  
453 mean OD<sub>450</sub> value at 500 pg/mL IL-6 are considered to be high responders. Low and high  
454 responders should be excluded from analysis. Each test must be performed with PBMCs  
455 from four different donors. From a set of four donors, a maximum of one donor may be  
456 excluded from the assessment; otherwise, the test must be repeated with four different  
457 donors.

## 458 9.0 DATA INTERPRETATION/DECISION CRITERIA

### 459 9.1 Decision Criteria for Determination of Pyrogenicity

460 The validity of the endotoxin standard curve should be calculated using a four-parameter  
461 logistic model. If necessary to satisfy the model, endotoxin concentrations may be modified.

462 Calculate the mean OD<sub>450</sub> values of all of the replicates in each experimental group. Calibrate  
463 the mean OD<sub>450</sub> value for each test substance using the endotoxin standard curve and  
464 document the estimated endotoxin concentration. Multiply the estimated endotoxin  
465 concentration by the dilution factor, if necessary. This value represents the pyrogenicity of  
466 the sample in terms of endotoxin equivalents for that particular donor. The *t*-test is used to  
467 compare the data of a test sample against the data of the EC (0.5 EU/mL) that is performed in  
468 parallel. If this test results in a significant *p*-value (i.e., smaller than 1%), then the sample is  
469 considered to be non-pyrogenic, and as pyrogenic if otherwise (Hoffmann et al., 2005), as  
470 long as the assay acceptability criteria in **Section 8.0** has been met.

#### 471 9.1.1 Decision Level 1

- 472 • If all donors show a negative reaction, then the product passes.
- 473 • If two or more donors show a positive reaction, then the product fails.
- 474 • If one donor shows a positive reaction, then an additional test with four donors  
475 has to be performed (go to Decision Level 2).

#### 476 9.1.2 Decision Level 2

- 477 • If out of 6 to 8 donors, only one donor shows a positive reaction, then the  
478 product passes.
- 479 • In any other case, the product fails.

## 480 10.0 STUDY REPORT

481 The test report should include the following information:

### 482 *Test Substances and Control Substances*

- 483 • Name of test substance



- 484 • Purity and composition of the substance or preparation
- 485 • Physicochemical properties (e.g., physical state, water solubility)
- 486 • Treatment of the test/control substances prior to testing (e.g., vortexing,
- 487 sonication, warming, resuspension solvent)

488 *Justification of the In Vitro Test Method and Protocol Used*

489 *Test Method Integrity*

- 490 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
- 491 test method over time
- 492 • If the test method employs proprietary components, documentation on the
- 493 procedure used to ensure their integrity from “lot-to-lot” and over time
- 494 • The procedures that the user may employ to verify the integrity of the
- 495 proprietary components

496 *Criteria for an Acceptable Test*

- 497 • Acceptable concurrent positive control ranges based on historical data
- 498 • Acceptable negative control data

499 *Test Conditions*

- 500 • Cell system used
- 501 • Calibration information for the spectrophotometer used to read the ELISA
- 502 • Details of test procedure used
- 503 • Description of any modifications of the test procedure
- 504 • Reference to historical data of the model
- 505 • Description of evaluation criteria used

506 *Results*

- 507 • Tabulation of data from individual test samples

508 *Description of Other Effects Observed*

509 *Discussion of the Results*

510 *Conclusion*

511 *A Quality Assurance Statement for Good Laboratory Practice (GLP)-Compliant Studies*

- 512           • This statement should indicate all inspections made during the study and the  
513           dates any results were reported to the Study Director. This statement should  
514           also confirm that the final report reflects the raw data.

515 If GLP-compliant studies are performed, then additional reporting requirements provided in  
516 the relevant guidelines (e.g., OECD 1998; EPA 2003a, 2003b; FDA 2003) should be  
517 followed.

518

518 **11.0 REFERENCES**

519 Bleeker WK, de Groot EM, den Boer PJ, Biessels PT, Aarden LA, Bakker JC. 1994.

520 Measurement of interleukin-6 production by monocytes for *in vitro* safety testing of

521 hemoglobin solutions. *Artif Cells Blood Substit Immobil Biotechnol* 22:835-840.

522 Dinarello CA, O'Connor JV, LoPreste G, Swift RL. 1984. Human leukocyte pyrogen test for

523 detection of pyrogenic material in growth hormone produced by recombinant *Escherichia*

524 *coli*. *J Clin Microbiol* 20:323-329.

525 EPA. 2003a. Good Laboratory Practice Standards. Toxic Substances Control Act. 40 CFR

526 792.

527 EPA. 2003b. Good Laboratory Practice Standards. Federal Insecticide, Fungicide, and

528 Rodenticide Act. 40 CFR 160.

529 FDA. 2003. Good Laboratory Practices for Nonclinical Laboratory Studies. 21 CFR 58.

530 Hoffmann S, Peterbauer A, Schindler S, Fennrich S, Poole S, Mistry Y, Montag-Lessing T,

531 Spreitzer I, Löschner B, van Aalderen M, Bos R, Gommer M, Nibbeling R, Werner-

532 Felmayer G, Loitzl P, Jungi T, Brcic M, Brügger P, Frey E, Bowe G, Casado J, Coecke S, de

533 Lange J, Mogster B, Næss LM, Aaberge IS, Wendel A, Hartung T. 2005. International

534 validation of novel pyrogen tests based on human monocytoïd cells. *J Immunol Methods*

535 298:161-173.

536 Organization for Economic Co-operation and Development (OECD). 1998. OECD Series on

537 Principle of Good Laboratory Practice and Compliance Monitoring. No. 1. OECD Principles

538 of Good Laboratory Practice (as revised in 1997). OECD, ENV/MC/CHEM (98)17. Paris:

539 OECD.

540 Poole S, Mistry Y, Ball C, Gaines Das RE, Opie LP, Tucker G, Patel M. 2003. A rapid 'one-

541 plate' *in vitro* test for pyrogens. *J Immunol Methods* 274:209-220.

542

542 **12.0 TERMINOLOGY AND FORMULA**

543 **12.1 Assay Sensitivity ( $\lambda$ )**

544 For an *in vitro* cell-based assay, the variable  $\lambda$  is defined as the lowest statistically significant  
545 point on the standard endotoxin concentration-response curve and represents the relative  
546 sensitivity of the test method for the detection of endotoxin (i.e., level of detection).

547 **12.2 Endotoxin Control (EC)**

548 The EC is incubated with PBMCs and serves as the positive control for the experiment. The  
549 results should be compared to historical values to insure that it provides a known level of  
550 cytokine release relative to the NSC.

551 **12.3 Endotoxin Limit Concentration (ELC)**

552 The ELC is the maximum allowable concentration of endotoxin for a particular product and  
553 is expressed in EU per volume (mL) or weight (mg). The ELC is defined by the FDA or  
554 specified in the USP<sup>2</sup>. It is calculated as the product of K/M, where:

555 K is the threshold pyrogen dose for parenteral use in rabbits or humans (5.0 EU/kg). At an  
556 injection volume of 10 mL/kg, K is equal to 0.5 EU/mL.

557 M is the larger of the rabbit dose or the maximum human dose administered in one hour as  
558 defined below and varies with test substance<sup>3</sup>.

559 **12.4 Maximum Valid Dilution (MVD)**

560 The MVD is the maximum dilution of a test substance that can be tolerated in a test system  
561 without exceeding the ELC, if the test substance must be diluted as a result of assay  
562 interference. Dilutions beyond the MVD would not be valid for endotoxin detection in the

---

<sup>2</sup> ELC values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, Pharmacopoeial guidelines, FDA publications).

<sup>3</sup> Values for most marketed pharmaceutical products are provided in the USP or in other pharmacopoeial or regulatory publications (e.g., European Pharmacopoeia, Japanese Pharmacopoeia, pharmacopoeial guidelines, FDA publications).

563 test system. Calculation of the MVD is dependent on whether or not the ELC for a test  
564 substance is published. When the ELC is known, the MVD is:

565 
$$\text{MVD} = (\text{ELC} \times \text{Product Potency [PP]})/\lambda$$

566 As an example, for "Cyclophosphamide Injection," the ELC is 0.17 EU/mg, PP is 20  
567 mg/mL, and the assay sensitivity is 0.1 EU/mL. The calculated MVD would be 34. The test  
568 substance can be diluted no more than 1:34 prior to testing.

569 If the ELC is not known, the MVD is:

570 
$$\text{MVD} = \text{PP}/\text{Minimum Valid Concentration (MVC)}$$

571 where,  $\text{MVC} = (\lambda \times \text{M})/\text{K}$

572 where, M is the maximum human dose

573 As an example, for "Cylophosphamide Injection," the PP is 20 mg/mL, M is 30  
574 mg/kg, and assay sensitivity is 0.1 EU/mL. The calculated MVC is 0.6 mg/mL and the MVD  
575 is 33.3. The test substance can be diluted no more than 1:33 in the assay prior to testing.

576 **12.5 Negative Product Control (NPC)**

577 The NPC is a test sample to which PFS is added. The NPC is the baseline for determination  
578 of cytokine release relative to the endotoxin-spiked PPC.

579 **12.6 Negative Saline Control (NSC)**

580 The NSC is PBMCs (in RPMI-C) incubated with PFS (used for dilution of test substance)  
581 and is used as the blank.

582 **12.7 Parenteral Threshold Pyrogen Dose (K)**

583 The value K represents the threshold pyrogen dose for parenteral products for rabbits and  
584 humans. Based on experimental data, K is fixed at 5.0 EU/kg. For intrathecal products, K is  
585 0.2 EU/kg.

586 **12.8 Positive Product Control (PPC)**

587 The PPC is a test substance spiked with the control standard endotoxin (i.e., 0.5 EU/mL or an  
588 amount of endotoxin equal to that which produces  $\frac{1}{2}$  the maximal increase in OD from the

589 endotoxin standard curve) to insure that the test system is capable of endotoxin detection in  
590 the product as diluted in the assay.

591 **12.9 Product Potency (PP)**

592 The concentration for a test substance is the PP typically expressed as  $\mu\text{g/mL}$  or  $\text{mg/mL}$ .

593 **12.10 Rabbit Pyrogen Test (RPT) Dose or Maximum Human Dose (M)**

594 The variable M represents the rabbit test dose or the maximum human dose in 1 hr. The  
595 variable M is expressed in  $\text{mg/kg}$  and varies with the test substance. For  
596 radiopharmaceuticals, M should be adjusted to account for product activity (radioactive  
597 decay) at time administration. An average human standard weight of 70 kg is used for the  
598 calculation. If a pediatric dose should be used and it is higher than the adult dose, then it  
599 should be used in the formula.

600

601